

1 **Supplemental Material**

2
3 ANONYMOUS AUTHOR(S)

4
5 Additional Key Words and Phrases: Chemistry, Domain Specific Language, Type system, Compiler, Microflu-
6 idics, Laboratory-on-a-Chip (LoC), Elecrowetting-on-Dielectric (EWoD)
7

8 CONTENTS

9

10	Contents	1
11	1 Typing Judgments	3
12	1.1 Syntax	3
13	1.2 Evaluation Rules	4
14	1.3 Type Checking Rules	6
15	1.4 Type Inference Rules	8
16	2 Proofs	10
17	2.1 Helper Lemmas	10
18	2.2 Proof of Progress	14
19	2.3 Proof of Preservation	18
20	2.4 Proof of Soundness	24
21	2.5 Proof of Completeness	27
22	3 Syntax Study: ELISA Protocols	31
23	4 Background	33
24	4.1 Digital Microfluidic Technology	33
25	5 Assay Execution Videos	35
26	6 BioScript’s ChemType Tests	35
27	6.1 Real-world Prevention Tests	35
28	6.2 Synthetic Failures	36
29	6.3 Synthetic Passes	37
30	7 BioScript Assays	37
31	7.1 Syntactic Sugar Translation	37
32	7.2 Urine Opiate Hierarchy	37
33	7.3 PCR Droplet Replenishment	41
34	7.4 Probabilistic PCR	41
35	7.5 ChemType Test1	42
36	7.6 ChemType Test2	42
37	7.7 ChemType Test3	43
38	7.8 Broad Spectrum Opiate	43
39	7.9 Cipro ELISA	44
40	7.10 Diazepam ELISA	44
41	7.11 547 Dilution	45
42	7.12 Fentanyl ELISA	45
43	7.13 Morphine ELISA	46
44	7.14 Morphine ELISA with Control Samples	47
45	7.15 Heroin ELISA	48
46	7.16 Oxycodone ELISA	48

50	7.17	Image Probe Synthesis	49
51	7.18	Glucose Detection	49
52	7.19	PCR	49
53	7.20	Neurotransmitter Sensing	50
54	8	AquaCore Assays	50
55	8.1	Glucose Detection	50
56	8.2	PCR	51
57	8.3	Imaging Probe Synthesis	51
58	8.4	Neurotransmitter Sensing	52
59	9	BioCoder Assays	53
60	9.1	PCR	53
61	9.2	Probabilistic PCR	53
62	9.3	PCR Droplet Replenish	55
63	9.4	Glucose Detection	56
64	9.5	Image Probe Synthesis	57
65	9.6	Neurotransmitter Sensing	58
66	10	Antha Assays	59
67	10.1	Glucose Detection	59
68	10.2	Imaging Probe Synthesis	61
69	10.3	PCR	62
70	10.4	Neurotransmitter Sensing	66
71	11	ChemType Results	68
72	11.1	Efficacy of ChemType	68
73	11.2	ChemType Runtime Evaluation	68
74		References	69
75			
76			
77			
78			
79			
80			
81			
82			
83			
84			
85			
86			
87			
88			
89			
90			
91			
92			
93			
94			
95			
96			
97			
98			

1 TYPING JUDGMENTS

Typing judgments are repeated here for convenience

1.1 Syntax

$t ::=$		Terms:
	$x \in \mathcal{X}$	Variable
	$t_1 \oplus t_2$	Math operation
	detect module on x for t	Detect
	$v \in \mathcal{V}$	Value
$v ::=$		Values:
	mat	Material value
	r	Real number
	n	Natural number
$module ::=$		
	$module_1, \dots, module_n$	Sensor module(s)
$s ::=$		Statements:
	$i; s$	Sequencing
	skip	Skip
$i ::=$		instructions:
	$x := t$	Assignment
	$x := \text{mix } x_1 \text{ with } x_2 \text{ for } t$	Mixing
	$\langle x_1, \dots, x_n \rangle := \text{split } x \text{ into } n$	Splitting
	if t then s_1 else s_2	Conditional
	while t s	Loop
$T ::=$		Union Types:
	$\cup \bar{S}$	Union type
	V	Type variables
$S ::=$		Scalar types:
	$Mat_1 \mid \dots \mid Mat_n$	Material types
	\mathbb{R}	Real number
	\mathbb{N}	Natural number
$\Gamma ::=$		Context:
	\emptyset	Empty context
	$\Gamma, x : T$	Variable type binding
\mathcal{X}		Set of variables x
\mathcal{C}		Constraints

Fig. 1. Syntax of ChemType.

1.2 Evaluation Rules

$$\begin{array}{l}
 \text{E-VAR} \\
 \frac{x \in \text{dom}(\sigma)}{(\sigma, x) \rightarrow \sigma(x)} \\
 \\
 \text{E-MATHR1} \\
 \frac{(\sigma, t_1) \rightarrow t'_1}{(\sigma, t_1 \oplus t_2) \rightarrow t'_1 \oplus t_2} \\
 \\
 \text{E-MATHR2} \\
 \frac{v \in \mathbb{N} \vee v \in \mathbb{R} \quad (\sigma, t_2) \rightarrow t'_2}{(\sigma, v \oplus t_2) \rightarrow v \oplus t'_2} \\
 \\
 \text{E-MATH} \\
 \frac{(v_1 \in \mathbb{N} \wedge v_2 \in \mathbb{N}) \vee (v_1 \in \mathbb{R} \wedge v_2 \in \mathbb{R}) \quad v_1 \oplus v_2 = v}{(\sigma, v_1 \oplus v_2) \rightarrow v} \\
 \\
 \text{E-DETECTR} \\
 \frac{(\sigma, t) \rightarrow t'}{(\sigma, \text{detect } module \text{ on } x \text{ for } t) \rightarrow \text{detect } module \text{ on } x \text{ for } t'} \\
 \\
 \text{E-DETECT} \\
 \frac{\sigma(x) \in Mat_i \quad r_2 = \text{detect}(\sigma(x), module, r_1)}{(\sigma, \text{detect } module \text{ on } x \text{ for } r_1) \rightarrow r_2}
 \end{array}$$

Fig. 2. Evaluation rules for terms.

197	E-ASSIGNR		E-ASSIGN'
198	$(\sigma, t) \rightarrow t' \quad t \notin \mathcal{X}$	$\frac{\text{E-ASSIGN}}{(\sigma, x := v; s) \rightarrow (\sigma[x \mapsto v], s)}$	$\frac{\sigma' = (\sigma \setminus \{x'\})[x \mapsto \sigma(x')]}{(\sigma, x := x'; s) \rightarrow (\sigma', s)}$
199	$(\sigma, x := t; s) \rightarrow (\sigma, x := t'; s)$		
200			
201	E-MixR		
202	$(\sigma, t) \rightarrow t'$		
203	$(\sigma, x := \text{mix } x_1 \text{ with } x_2 \text{ for } t; s) \rightarrow (\sigma, x := \text{mix } x_1 \text{ with } x_2 \text{ for } t'; s)$		
204			
205	E-Mix		
206	$\sigma(x_1) \in \text{Mat}_i \quad \sigma(x_2) \in \text{Mat}_j$		
207	$\text{interact}(\sigma(x_1), \sigma(x_2), r) \neq \perp$		
208	$\sigma' = (\sigma \setminus \{x_1, x_2\})[x \mapsto \text{interact}(\sigma(x_1), \sigma(x_2), r)]$		
209	$(\sigma, x := \text{mix } x_1 \text{ with } x_2 \text{ for } r; s) \rightarrow (\sigma', s)$		
210			
211	E-SPLIT		
212	$\sigma(x) \in \text{Mat}_i \quad \sigma' = (\sigma \setminus \{x\})[x_i \mapsto \text{split}(\sigma(x), n)]$		
213	$(\sigma, \langle x_0, \dots, x_n \rangle := \text{split } x \text{ into } n; s) \rightarrow (\sigma', s)$		
214			
215			
216	E-IFR		
217	$(\sigma, t) \rightarrow t'$		
218	$(\sigma, \text{if } t \text{ then } s_1 \text{ else } s_2; s) \rightarrow (\sigma, \text{if } t' \text{ then } s_1 \text{ else } s_2; s)$		
219			
220	E-IFTRUE		
221	$n \neq 0$		
222	$(\sigma, \text{if } n \text{ } s_1 \text{ else } s_2; s) \rightarrow (\sigma, s_1 \bullet s)$		
223			
224	E-IFFALSE		
225	$(\sigma, \text{if } 0 \text{ } s_1 \text{ else } s_2; s) \rightarrow (\sigma, s_2 \bullet s)$		
226			
227	E-WHILE		
228	$(\sigma, \text{while } t \text{ } s_1; s_2) \rightarrow (\sigma, \text{if } t \text{ then } (s_1 \bullet \text{while } t \text{ } s_1; s_2) \text{ else } s_2)$		
229			
230	$\text{skip} \bullet s = s$		
231	$(i; s) \bullet s' = i; (s \bullet s')$		
232			

Fig. 3. Evaluation rules for statements.

1.3 Type Checking Rules

(T-VAR)	$\frac{x : T \in \Gamma \quad x \in X}{\Gamma, X \vdash x : T}$
(T-MATH)	$\frac{\Gamma, X \vdash t_1 : T \quad \Gamma, X \vdash t_2 : T \quad T = \mathbb{N} \vee T = \mathbb{R}}{\Gamma, X \vdash t_1 \oplus t_2 : T}$
(T-DETECT)	$\frac{\Gamma, X \vdash x : \cup \overline{Mat_i} \quad \Gamma, X \vdash t : \mathbb{R}}{\Gamma, X \vdash \text{detect module on } x \text{ for } t : \mathbb{R}}$
(T-MAT)	$\frac{\overline{mat \in Mat_i}}{\Gamma, X \vdash mat : \cup \overline{Mat_i}}$
(T-REAL)	$\Gamma, X \vdash r : \mathbb{R}$
(T-NAT)	$\Gamma, X \vdash n : \mathbb{N}$

Fig. 4. Type checking rules for terms.

(T-INST)	$\frac{\Gamma, X \vdash i, X' \quad \Gamma, X' \vdash s, X''}{\Gamma, X \vdash i; s, X''}$
(T-SKIP)	$\frac{}{\Gamma, X \vdash \text{skip}, X}$
(T-ASSIGN-1)	$\frac{x : T \in \Gamma \quad \Gamma, X \vdash v : T' \quad T' \subseteq T}{\Gamma, X \vdash x := v, X \cup \{x\}}$
(T-ASSIGN-2)	$\frac{x : T \in \Gamma \quad \Gamma, X \vdash x' : T' \quad T' \subseteq T}{\Gamma, X \vdash x := x', X \setminus \{x'\} \cup \{x\}}$
(T-ASSIGN-3)	$\frac{x : T \in \Gamma \quad t \notin \mathcal{V} \cup \mathcal{X} \quad \Gamma, X \vdash t : T' \quad T' = \mathbb{R} \vee T' = \mathbb{N} \quad T' \subseteq T}{\Gamma, X \vdash x := t, X \cup \{x\}}$
(T-MIX)	$\frac{\Gamma, X \vdash x_1 : \cup \overline{Mat_i} \quad \Gamma, X \vdash x_2 : \cup \overline{Mat_j} \quad \Gamma, X \vdash t : \mathbb{R} \quad \text{interact-abs}(Mat_i, Mat_j) \subseteq \Gamma(x)}{\Gamma, X \vdash x := \text{mix } x_1 \text{ with } x_2 \text{ for } t, X \setminus \{x_1, x_2\} \cup \{x\}}$
(T-SPLIT)	$\frac{\Gamma, X \vdash x : \cup \overline{Mat_i} \quad \Gamma(x) \subseteq \Gamma(x_1), \dots, \Gamma(x) \subseteq \Gamma(x_n)}{\Gamma, X \vdash \langle x_1, \dots, x_n \rangle := \text{split } x \text{ into } n, X \setminus \{x\} \cup \{x_1, \dots, x_n\}}$
(T-IF)	$\frac{\Gamma, X \vdash t : \mathbb{N} \quad \Gamma, X \vdash s_1, X' \quad \Gamma, X \vdash s_2, X''}{\Gamma, X \vdash \text{if } t \text{ then } s_1 \text{ else } s_2, X' \cap X''}$
(T-WHILE)	$\frac{\Gamma, X \vdash t : \mathbb{N} \quad \Gamma, X \vdash s, X' \quad X \subseteq X'}{\Gamma, X \vdash \text{while } t \text{ s}, X}$

Fig. 5. Type checking rules for statements.

1.4 Type Inference Rules

(CT-VAR)	$\frac{x : T \in \Gamma \quad x \in X}{\Gamma, X \vdash x : T \mid \emptyset}$
(CT-MATH)	$\frac{\Gamma, X \vdash t_1 : T_1 \mid C_1 \quad \Gamma, X \vdash t_2 : T_2 \mid C_2}{\Gamma, X \vdash t_1 \oplus t_2 : T_1 \mid C_1 \cup C_2 \cup \{T_1 = T_2 = \mathbb{N} \vee T_1 = T_2 = \mathbb{R}\}}$
(CT-DETECT)	$\frac{\Gamma, X \vdash x : T_1 \mid C_1 \quad \Gamma, X \vdash t : T_2 \mid C_2}{\Gamma, X \vdash \text{detect module on } x \text{ for } t : \mathbb{R} \mid C_1 \cup C_2 \cup \{T_1 \cap \{\mathbb{R}, \mathbb{N}\} = \emptyset, T_2 = \mathbb{R}\}}$
(CT-MAT)	$\frac{\overline{mat \in Mat_i}}{\Gamma, X \vdash mat : \cup \overline{Mat_i} \mid \emptyset}$
(CT-REAL)	$\Gamma, X \vdash r : \mathbb{R} \mid \emptyset$
(CT-NAT)	$\Gamma, X \vdash n : \mathbb{N} \mid \emptyset$

Fig. 6. Type Inference rules for terms.

(CT-INST)	$\frac{\Gamma, X \vdash i, X' \mid C_1 \quad \Gamma, X' \vdash s, X'' \mid C_2}{\Gamma, X \vdash i; s, X'' \mid C_1 \cup C_2}$
(CT-SKIP)	$\frac{}{\Gamma, X \vdash \text{skip}, X \mid \emptyset}$
(CT-ASSGN-1)	$\frac{x : T \in \Gamma \quad \Gamma, X \vdash v : T' \mid C'}{\Gamma, X \vdash x := v, X \cup \{x\} \mid C' \cup \{T' \subseteq T\}}$
(CT-ASSGN-2)	$\frac{x : T \in \Gamma \quad \Gamma, X \vdash x' : T' \mid C'}{\Gamma, X \vdash x := x', X \setminus \{x'\} \cup \{x\} \mid C' \cup \{T' \subseteq T\}}$
(CT-ASSGN-3)	$\frac{x : T \in \Gamma \quad t \notin \mathcal{V} \cup \mathcal{X} \quad \Gamma, X \vdash t : T' \mid C'}{\Gamma, X \vdash x := t, X \cup \{x\} \mid C' \cup \{T' = \mathbb{R} \vee T' = \mathbb{N}, T' \subseteq T\}}$
(CT-MIX)	$\frac{\Gamma, X \vdash x_1 : T \mid C \quad \Gamma, X \vdash x_2 : T' \mid C' \quad \Gamma, X \vdash t : T'' \mid C''}{\Gamma, X \vdash x := \text{mix } x_1 \text{ with } x_2 \text{ for } t, X \setminus \{x_1, x_2\} \cup \{x\} \mid C \cup C' \cup C'' \cup \{T \cap \{\mathbb{R}, \mathbb{N}\} = \emptyset, T' \cap \{\mathbb{R}, \mathbb{N}\} = \emptyset, T'' = \mathbb{R}, \text{Mat}_i \in T \wedge \text{Mat}_j \in T' \Rightarrow \text{interact-abs}(\text{Mat}_i, \text{Mat}_j) \subseteq \Gamma(x)\}}$
(CT-SPLIT)	$\frac{\Gamma, X \vdash x : T \mid C}{\Gamma, X \vdash \langle x_1, \dots, x_n \rangle := \text{split } x \text{ into } n, X \setminus \{x\} \cup \{x_1, \dots, x_n\} \mid C \cup \{T \cap \{\mathbb{R}, \mathbb{N}\} = \emptyset, T \subseteq \Gamma(x_1), \dots, T \subseteq \Gamma(x_n)\}}$
(CT-IF)	$\frac{\Gamma, X \vdash t : T \mid C \quad \Gamma, X \vdash s_1, X' \mid C_1 \quad \Gamma, X \vdash s_2, X'' \mid C_2}{\Gamma, X \vdash \text{if } t \text{ then } s_1 \text{ else } s_2, X' \cap X'' \mid C \cup C_1 \cup C_2 \cup \{T = \mathbb{N}\}}$
(CT-WHILE)	$\frac{\Gamma, X \vdash t : T \mid C \quad \Gamma, X \vdash s, X' \mid C' \quad X \subseteq X'}{\Gamma, X \vdash \text{while } t \text{ s}, X, C \cup C' \cup \{T = \mathbb{N}\}}$

Fig. 7. Type Inference rules for statements.

2 PROOFS

2.1 Helper Lemmas

LEMMA 1 (NAME EXTENSION).

$\forall \Gamma, X, i, X', X''.$

$\Gamma, X_1 \vdash i, X_2 \wedge X_1 \subseteq X'_1 \Rightarrow$

$\exists X'_2.$

$\Gamma, X'_1 \vdash i, X'_2 \wedge X_2 \subseteq X'_2$

and

$\forall \Gamma, X, s, X', X''.$

$\Gamma, X_1 \vdash s, X_2 \wedge X_1 \subseteq X'_1 \Rightarrow$

$\exists X'_2.$

$\Gamma, X'_1 \vdash s, X'_2 \wedge X_2 \subseteq X'_2$

PROOF. Trivial by mutual induction on s and i . □

LEMMA 2.

For every Γ, X, i, X', s , and X'' , if

$\Gamma, X \vdash i, X', X'' \subseteq X'$ and $\Gamma, X'' \vdash s, X'''$

there exists X'''' such that

$\Gamma, X \vdash i; s, X''''$ and $X''' \subseteq X''''$

PROOF. Direct from Lemma 1 and the rule T-INST. □

LEMMA 3 (TYPING \bullet).

For every Γ, X, s_1, X', s_2 , and X''

$\Gamma, X \vdash s_1, X' \wedge \Gamma, X' \vdash s_2, X'' \Rightarrow \Gamma, X \vdash s_1 \bullet s_2, X''$

PROOF. Trivial by induction on s_1 . □

LEMMA 4.

For every Γ, X, s_1, X', s_2 , and X'' , if

$\Gamma, X \vdash s_1, X', X'' \subseteq X'$ and $\Gamma, X'' \vdash s_2, X'''$

there exists X'''' such that

$\Gamma, X \vdash s_1 \bullet s_2, X''''$ and $X''' \subseteq X''''$

PROOF. Direct from Lemma 1 and Lemma 3. □

LEMMA 5 (CANONICAL FORMS).

For every Γ, X and v ,

- If $\Gamma, X \vdash v : \text{Mat}_i$, then $v \in \text{Mat}_i$.
- If $\Gamma, X \vdash v : \mathbb{R}$, then $v \in \mathbb{R}$.
- If $\Gamma, X \vdash v : \mathbb{N}$, then $v \in \mathbb{N}$.

PROOF. Immediate from case analysis on the structure of v and using the inversion lemma, Lemma 6. □

LEMMA 6 (INVERSION ON TYPING OF TERMS).

(1) If $\Gamma, X \vdash x : T$,

then $x : T \in \Gamma$ and $x \in X$.

(2) If $\Gamma, X \vdash t_1 \oplus t_2 : T$,

then $\Gamma, X \vdash t_1 : T$ and $\Gamma, X \vdash t_2 : T$ and $T = \mathbb{N} \vee T = \mathbb{R}$

- (3) If $\Gamma, X \vdash \text{detect module on } x \text{ for } t : T$,
 then there exists $\overline{\text{Mat}}_i$ such that $\Gamma, X \vdash x : \cup \overline{\text{Mat}}_i$ and $\Gamma, X \vdash t : \mathbb{R}$ and $T = \mathbb{R}$
- (4) If $\Gamma, X \vdash \text{mat} : T$, then there exists i such that $T = \text{Mat}_i$
- (5) If $\Gamma, X \vdash r : T$, then $T = \mathbb{R}$
- (6) If $\Gamma, X \vdash n : T$, then $T = \mathbb{N}$

PROOF. Immediate from case analysis on the type derivation rules. \square

LEMMA 7 (INVERSION ON TYPING OF STATEMENTS AND INSTRUCTIONS).

- (1) For all Γ, X, i, s, X'' ,
 If
 $\Gamma, X \vdash i; s, X''$,
 then, there exists X' such that
 $\Gamma, X \vdash i, X'$ and
 $\Gamma, X' \vdash s, X''$.
- (2) For all Γ, X, x, t, X' ,
 If
 $\Gamma, X \vdash x := t, X'$
 then, there exists T, T' such that
 $x : T \in \Gamma$,
 $\Gamma, X \vdash t : T'$,
 $T' \subseteq T$,
 $T' = \mathbb{R} \vee T' = \mathbb{N} \vee t = \text{mat}$, and
 $X' = X \cup \{x\}$
 or
 $t = x'$ and
 $X' = X \setminus \{x'\} \cup \{x\}$.
- (3) For all $\Gamma, X, x, x_1, x_2, t, X'$,
 If
 $\Gamma, X \vdash x := \text{mix } x_1 \text{ with } x_2 \text{ for } t, X'$
 then, there exist i and j such that
 $\Gamma, X \vdash x_1 : \cup \overline{\text{Mat}}_i$,
 $\Gamma, X \vdash x_2 : \cup \overline{\text{Mat}}_j$,
 $\Gamma, X \vdash t : \mathbb{R}$,
 $\text{interact-abs}(\text{Mat}_i, \text{Mat}_j) \subseteq \Gamma(x)$, and
 $X' = X \setminus \{x_1, x_2\} \cup \{x\}$.
- (4) For all $\Gamma, X, x, x_1, \dots, x_n, X'$,
 If
 $\Gamma, X \vdash \langle x_1, \dots, x_n \rangle := \text{split } x \text{ into } n, X'$
 then, there exist i such that
 $\Gamma, X \vdash x : \cup \overline{\text{Mat}}_i$,
 $\Gamma(x) \subseteq \Gamma(x_1), \dots, \Gamma(x) \subseteq \Gamma(x_n)$, and
 $X' = X \setminus \{x\} \cup \{x_1, \dots, x_n\}$.
- (5) For all $\Gamma, X, t, s_1, s_2, X'''$,
 If
 $\Gamma, X \vdash \text{if } t \text{ then } s_1 \text{ else } s_2, X'''$
 then, there exists X' and X'' such that
 $\Gamma, X \vdash t : \mathbb{N}$,

540 $\Gamma, X \vdash s_1, X'$,
 541 $\Gamma, X \vdash s_2, X''$, and
 542 $X''' = X' \cap X''$.
 543 (6) For all Γ, X, t, s, X' ,
 544 If
 545 $\Gamma, X \vdash \text{while } t \text{ } s, X'$
 546 then,
 547 $\Gamma, X \vdash t : \mathbb{N}$
 548 $\Gamma, X \vdash s, X''$,
 549 $X \subseteq X''$ and
 550 $X' = X$.

551 PROOF. Immediate from case analysis on the type derivation rules. □

552 LEMMA 8.

553 For every σ, t and t' if $(\sigma, t) \rightarrow t'$ then $t' \notin X$

554 PROOF. Immediate from the term transition rules. □

Helper Definitions

The conservative property of the abstract interact-abs function:

$\forall mat_i, mat_j, r.$

$mat_i \in Mat_i \wedge mat_j \in Mat_j \Rightarrow$

$interact(mat_i, mat_j, r) = \perp \Rightarrow interact-abs(Mat_i, Mat_j) \text{ undefined}$

$interact(mat_i, mat_j, r) \neq \perp \Rightarrow interact(mat_i, mat_j, r) \in interact-abs(Mat_i, Mat_j)$

The type of the functions used for evaluation:

$detect : Mat_i \rightarrow Module \rightarrow \mathbb{R} \rightarrow \mathbb{R}$

$interact : Mat_i \rightarrow Mat_j \rightarrow interact-abs(Mat_i, Mat_j)$

$split : Mat_i \rightarrow \mathbb{N} \rightarrow Mat_i$

We define the consistency condition between the static typing environment Γ and the runtime store σ as:

$consistent(\Gamma, X, \sigma) = \forall x, T.$

$(x : T) \in \Gamma \wedge x \in X \Rightarrow$

$\sigma(x) \in T$

2.2 Proof of Progress

LEMMA 9 (**PROGRESS OF TERMS**).

For every Γ, X, t or T ,

if

$\Gamma, X \vdash t : T$

then

$\forall \sigma. \text{consistent}(\Gamma, X, \sigma) \Rightarrow \exists t'. (\sigma, t) \rightarrow t'$ or
 t is a value.

PROOF.

Proof by induction on t and case analysis thereafter.

We assume

(1) $\Gamma, X \vdash t : T$

Case for the rule T-VAR:

We have

(2) $t = x$

(3) $x : T \in \Gamma$

(4) $x \in X$

We assume

(5) $\text{consistent}(\Gamma, X, \sigma)$

From [3], [4] and [5]:

(5) $\sigma(x) \in T$

By the rule E-VAR on [5]:

(6) $(\sigma, x) \rightarrow \sigma(x)$

By the rule E-VAR on [6] and [2]:

(6) $(\sigma, t) \rightarrow \sigma(x)$

Case for the rule T-MATH:

(2) $t = t_1 \oplus t_2$

(3) $\Gamma \vdash t_1 : T$

(4) $\Gamma \vdash t_2 : T$

(5) $T = \mathbb{R} \vee T = \mathbb{N}$

By I.H. on t_1 :

Case 1: t_1 is not a value:

(6) $\forall \sigma, X. \text{consistent}(\sigma, X, \Gamma) \Rightarrow \exists t'_1. (\sigma, t_1) \rightarrow t'_1$

We assume that

(7) $\text{consistent}(\sigma, X, \Gamma)$

and prove that

$\exists t', (\sigma, t) \rightarrow t'$

From [6] and [7], there exists t'_1 such that

(8) $(\sigma, t_1) \rightarrow t'_1$

From the rule E-MATHR1 on [8]:

(9) $(\sigma, t_1 \oplus t_2) \rightarrow t'_1 \oplus t_2$

From [9] on [2]

$(\sigma, t) \rightarrow t'_1 \oplus t_2$

Case 2:

(10) t_1 is a value v_1 :

By Lemma 5 on [3], [5], [10]

(11) $v_1 \in \mathbb{N} \vee v_1 \in \mathbb{R}$

By I.H. on t_2 :

Case 2.1: t_2 is not a value:

(12) $\forall \sigma, X. \text{consistent}(\sigma, X, \Gamma) \Rightarrow \exists t'_2, (\sigma, t_2) \rightarrow t'_2$

We assume that

(13) $\text{consistent}(\sigma, X, \Gamma)$

and prove that

$\exists t', (\sigma, t) \rightarrow t'$

From [12] and [13], there exists t'_2 such that

(14) $(\sigma, t_2) \rightarrow t'_2$

From the rule E-MATHR2 on [14]:

(15) $(\sigma, v_1 \oplus t_2) \rightarrow v_1 \oplus t'_2$

From [15], [2], [10]:

$(\sigma, t) \rightarrow v_1 \oplus t'_2$

Case 2.2:

(16) t_2 is a value, v_2 :

By Lemma 5 on [4], [5], [16]

(17) $v_2 \in \mathbb{N} \vee v_2 \in \mathbb{R}$

From the rule E-MATH on [11] and [17]:

(18) $(\sigma, v_1 \oplus v_2) \rightarrow v_1 \oplus v_2$

From on [17], [10], [16]

$(\sigma, t) \rightarrow v_1 \oplus v_2$

Case for the rule T-DETECT:

$t = \text{detect module}_i \text{ on } x \text{ for } t'$

Similar to the the rule T-MATH rule, by induction hypothesis in t' and then using the rule E-DETECT and the rule E-DETECTR. The consistency condition is used to derive the first premise of the rule E-DETECT.

Case for the rule T-MAT:

mat is a value.

Case for the rule T-REAL:

r is a value.

Case for the rule T-NAT:

n is a value. □

PROOF.

Case analysis on the typing derivation:

Case for the rule T-SKIP:

$s = \text{skip}$

Case for the rule T-INST:

$s = i; s'$

$\Gamma, X \vdash i, X''$

Immediate from Lemma 10. □

LEMMA 10 (**PROGRESS FOR INSTRUCTIONS**).

For every Γ, X, i, s, X' ,

if

$\Gamma, X \vdash i, X'$

then

$\forall \sigma. \text{consistent}(\Gamma, X, \sigma) \Rightarrow$

$\exists \sigma', s'. (\sigma, i; s) \rightarrow (\sigma', s')$

PROOF.

We have that

(1) $\Gamma, X \vdash i, X'$

Case analysis on the typing derivation:

Case for the rule T-ASSIGN-1:

(2) $i = (x := v)$

(3) $x : T \in \Gamma$

(4) $\Gamma, X \vdash t : T'$

(5) $T' \subseteq T$

From the rule E-ASSIGN

(6) $(\sigma, x := v; s) \rightarrow (\sigma[x \mapsto v]; s)$

From [6], [2]

$(\sigma, i; s) \rightarrow (\sigma[x \mapsto v]; s)$

Case for the rule T-ASSIGN-2:

Similar to the previous case. The reduction uses the rule E-ASSIGN'

Case for the rule T-ASSIGN-3:

(2) $i = (x := t)$

(3) $x : T \in \Gamma$

(4) $\Gamma, X \vdash t : T'$

(5) $(T' = \mathbb{R} \vee T' = \mathbb{N}) \wedge t \notin \mathcal{V} \cup \mathcal{X}$

(6) $T' \subseteq T$

By Lemma 9 on [4]:

Case 1:

(7) $\forall \sigma. \text{consistent}(\Gamma, X, \sigma) \Rightarrow \exists t'. (\sigma, t) \rightarrow t'$

We assume that

(8) $\text{consistent}(\Gamma, X, \sigma)$

We prove that

$(\sigma, i; s) \rightarrow (\sigma', s')$

From [7] and [8], there exists t' such that

(9) $(\sigma, t) \rightarrow t'$

From the rule E-ASSIGNR on [9] and [5]:

(10) $(\sigma, (x := t); s) \rightarrow (\sigma', (x := t'); s)$

From [2] on [10]:

$(\sigma, i; s) \rightarrow (\sigma', (x := t'); s)$

Case 2:

(11) t is a value, v .

Contradiction with [5].

Case for the rule T-Mix:

(2) $i = (x := \text{mix } x_1 \text{ with } x_2 \text{ for } t)$

The proof is similar to the case for the rule T-Assign. Lemma 9 is applied to the type derivation for t . There are two cases. Case 1: If t steps, the rule E-MixR is applicable. Case 2: If t is a value, the rule E-Mix is applicable. Lemma 6 and the consistency condition is used to show that x_1 and x_2 are both material values $\sigma(x_1)$ and $\sigma(x_2)$ in the store. In addition, from the case analysis on the typing derivation, we have that $\text{interact-abs}(Mat_i, Mat_j)$ is defined; thus, by the conservative property of the abstract interaction $\text{interact}(\sigma(x_1), \sigma(x_2)) \neq \perp$.

Case for the rule T-Split:

(2) $i = (\langle x_1, \dots, x_n \rangle = \text{split } x_1 \text{ into } n)$

The consistency condition is used to show that x is a material value in the store. Then, the rule E-Split is applicable.

Case for the rule T-If:

(2) $i = \text{if } t \text{ then } s_1 \text{ else } s_2$

We apply Lemma 9 to the type derivation for t . There are two cases. Case 1: If t steps, the rule E-IfR is applicable. Case 2: If t is a value, by Lemma 5 on the typing judgement for t , we know that it is a natural number. If it is non-zero, the rule E-IfTrue is applicable; otherwise the rule E-IfFalse is applicable.

Case for the rule T-While:

(2) $i = \text{while } t \text{ s}$

The rule E-While is applied without any premise.

□

LEMMA 11 (PROGRESS FOR STATEMENTS).

For every Γ, X, s, X' ,

if

$\Gamma, X \vdash s, X'$

then either s is skip or

$\forall \sigma. \text{consistent}(\Gamma, X, \sigma) \Rightarrow$

$\exists \sigma', s'. (\sigma, s) \rightarrow (\sigma', s')$

PROOF.

We have that

(1) $\Gamma, X \vdash s, X'$

Case analysis on s :

Case $s = \text{skip}$

Conclusion is immediate.

Case

(2) $s = i; s'$

From [1] and [2],
 (3) $\Gamma, X \vdash i; s', X'$
 By Lemma 7 on [3], there exists X'' such that
 (4) $\Gamma, X \vdash i, X''$
 (5) $\Gamma, X'' \vdash s', X'$
 By Lemma 10 on [4],
 (6) $\forall \sigma. \text{consistent}(\Gamma, X, \sigma) \Rightarrow$
 $\exists \sigma', s''. (\sigma, i; s') \rightarrow (\sigma', s'')$
 The conclusion is immediate from [2] and [6].

□

2.3 Proof of Preservation

LEMMA 12 (PRESERVATION OF TERMS).

For every Γ, X, t, T and σ ,

if

$\Gamma, X \vdash t : T$ and
 $(\sigma, t) \rightarrow t'$ and
 $\text{consistent}(\Gamma, X, \sigma)$

then

$\Gamma, X \vdash t' : T$

PROOF.

We have

(1) $\Gamma, X \vdash t : T$
 (2) $(\sigma, t) \rightarrow t'$
 (3) $\text{consistent}(\Gamma, X, \sigma)$

Straightforward induction on the derivation of $\Gamma, X \vdash t : T$ and then case analysis on the final rule in the derivation of $(\sigma, t) \rightarrow t'$

Case for the rule T-VAR:

From the rule T-VAR:

(4) $t = x$
 (5) $x : T \in \Gamma$
 (6) $x \in X$

From the rule E-VAR:

(7) $t' = \sigma(x)$

From [3], [5], and [6]

(7) $\sigma(x) \in T$

By case analysis on the value $\sigma(x)$ and the rule T-MAT,
 the rule T-REAL and the rule T-NAT

$\Gamma, X \vdash \sigma(x) : T$

Case for the rule T-MATH:

(4) $t = t_1 \oplus t_2$
 (5) $\Gamma, X \vdash t_1 : T$
 (6) $\Gamma, X \vdash t_2 : T$
 (7) $T = \mathbb{R} \vee T = \mathbb{N}$

Case analysis on [2]:

Case for the rule E-MATHR1:

$$(8) t' = t'_1 \oplus t_2$$

$$(9) (\sigma, t_1) \rightarrow t'_1$$

By I.H. on [5], [9], and [3]:

$$(10) \Gamma, X \vdash t'_1 : T$$

From the rule T-MATH on [8], [10], [6], [7]:

$$(11) \Gamma, X \vdash t' : T$$

Case for the rule E-MATHR2:

$$(12) t' = v_1 \oplus t'_2$$

$$(13) (\sigma, t_2) \rightarrow t'_2$$

By I.H. on [6], [13], and [3]:

$$(14) \Gamma, X \vdash t'_2 : T$$

From the rule T-MATH on [12], [5], [14]:

$$(15) \Gamma, X \vdash t' : T$$

Case for the rule E-MATH:

$$(16) t_1 = v_1$$

$$(17) t_2 = v_2$$

$$(18) (v_1 \in \mathbb{N} \wedge v_2 \in \mathbb{N}) \vee (v_1 \in \mathbb{R} \wedge v_2 \in \mathbb{R})$$

$$(19) v_1 \oplus v_2 = v$$

$$(20) t' = v$$

From [18], we consider the case:

$$(21) v_1 \in \mathbb{N} \wedge v_2 \in \mathbb{N}$$

The other case is similar.

From [21], [19]:

$$(22) v \in \mathbb{N}$$

From [20], [22] and the rule T-NAT:

$$(23) \Gamma, X \vdash t' : \mathbb{N}$$

From Lemma 6 on [5], [16] and [21]

$$(24) T = \mathbb{N}$$

From [23] and [24]:

$$(23) \Gamma, X \vdash t' : T$$

Case for the rule T-DETECT:

$$(4) t = \text{detect } module_i \text{ on } x \text{ for } t'$$

We consider the two cases for [2]: Case for the rule E-DETECTR: Induction hypothesis is applied to t' and then the rule T-DETECT is applied. Case for the rule E-DETECT: Immediate from the rule T-REAL.

Case for the rule T-MAT:

$$(4) t = mat$$

mat is a value and does not step.

Case for the rule T-REAL:

$$(4) t = r$$

r is a value and does not step.

Case for the rule T-NAT:

(4) $t = n$

n is a value and does not step.

□

LEMMA 13 (**PRESERVATION OF STATEMENTS**).

For every $\Gamma, X, \sigma, s, X'', \sigma', s'$,

if

$\Gamma, X \vdash s, X''$ and

$(\sigma, s) \rightarrow (\sigma', s')$ and

$\text{consistent}(\Gamma, X, \sigma)$

then there exists X' such that

$\Gamma, X' \vdash s', X''$ and

$\text{consistent}(\Gamma, X', \sigma')$

PROOF.

We have

(1) $\Gamma, X \vdash s, X''$

(2) $(\sigma, s) \rightarrow (\sigma', s')$

(3) $\text{consistent}(\Gamma, X, \sigma)$

Case Analysis on [1]:

Case for the rule T-SKIP:

Contradiction in [2]: The statement skip does not step.

Case for the rule T-INST:

(4) $\Gamma, X \vdash i, X'$

(5) $\Gamma, X' \vdash s'', X''$

(6) $s = i; s''$

Case Analysis on [4]:

Case for the rule T-ASSIGN-1:

(7) $i = (x := v)$

(8) $x : T \in \Gamma$

(9) $\Gamma, X \vdash t : T'$

(10) $X' = X \cup \{x\}$

(11) $T' \subseteq T$

From [6], [7]:

(12) $s = (x := v); s''$

Case analysis on [2]:

Case for the rule E-ASSIGN:

(13) $t = v$

(14) $s' = s''$

Case analysis on v :

Case (15) $v = r$

By Lemma 6 on [9], [13] and [15]

(16) $T' = \mathbb{R}$

From [16], [11]

(17) $\{\mathbb{R}\} \subseteq T$
 From [15], [17]
 (18) $v \in T$
 From [3], [10], [8], [18]
 (19) $\text{consistent}(\Gamma, X', \sigma[x \mapsto v])$
 From [5], [14]
 (20) $\Gamma, X' \vdash s', X''$
 The conclusion is [20] and [19].
 Case $T' = \mathbb{N}$
 Similar to the previous case.
 Case $t = \text{mat}$
 Similar to the previous case.
 Case for the rule E-ASSIGNR:
 $(\sigma, v) \rightarrow t'$
 There is no reduction rule for values. Contradiction.
 Case for the rule T-ASSIGN-2:
 The reduction is with the the rule E-ASSIGN'. Similar to the case for the rule T-ASSIGN-2, the consistency of σ with Γ and X and that $T' \subseteq T$ implies the consistency of $(\sigma \setminus \{x'\})[x \mapsto \sigma(x')]$ with Γ and $X \setminus \{x'\} \cup \{x\}$.
 Case for the rule T-ASSIGN-3:
 (7) $i = (x := t)$
 (8) $x : T \in \Gamma$
 (9) $\Gamma, X \vdash t : T'$
 (10) $(T' = \mathbb{R} \vee T' = \mathbb{N}) \wedge t \notin \mathcal{V} \cup \mathcal{X}$
 (11) $T' \subseteq T$
 (12) $X' = X \cup \{x\}$
 From [6], [7]:
 (13) $s = (x := t); s''$
 Case analysis on [2]:
 Case for the rule E-ASSIGNR:
 (14) $s' = (x := t'); s''$
 (15) $(\sigma, t) \rightarrow t'$
 By Lemma 12 on [9], [15], and [3]:
 (16) $\Gamma, X \vdash t' : T$
 By Lemma 8 on [15]
 (17) $t' \notin \mathcal{X}$
 From [17] and either the rule T-ASSIGN-1 or the rule T-ASSIGN-3 on [8], [16], [10], [11]:
 (18) $\Gamma, X \vdash x := t', X \cup \{x\}$
 From [18] and [12]:
 (19) $\Gamma, X \vdash x := t', X'$
 From the rule T-INST on [19], [5], and then [14]:
 (20) $\Gamma, X \vdash s', X''$
 The conclusion is [20], [3].
 Case for the rule E-ASSIGN:

$t = v$

Contradiction with [10].

Case for the rule E-Assign':

$t = x$

Contradiction with [10].

Case for the rule T-Mix:

(7) $i = (x := \text{mix } x_1 \text{ with } x_2 \text{ for } t)$

(8) $\Gamma, X \vdash x_1 : \cup \overline{Mat_i}$

(9) $\Gamma, X \vdash x_2 : \cup \overline{Mat_j}$

(10) $\Gamma, X \vdash t : \mathbb{R}$

(11) $\text{interact-abs}(Mat_i, Mat_j) \subseteq \Gamma(x)$

(12) $X \setminus \{x_1, x_2\} \cup \{x\}$

Case analysis on [2]: There are two cases.

Case of the rule E-MixR:

By Lemma 12, the type of t is preserved for t' . Thus, the rule T-Mix is applied to the new mix instruction. The assumed consistency condition is preserved since the store σ stays unchanged in the step.

Case of the rule E-Mix:

We already have the typing judgement for the remaining statement s'' in [5]. By Lemma 6 on [8] and the consistency condition, there exists i such that $\sigma(x_1) \in Mat_i$ where $\Gamma(x_1) = \cup \overline{Mat_i}$. Similarly, there exists j , $\sigma(x_2) \in Mat_j$ where $\Gamma(x_2) = \cup \overline{Mat_j}$. Since $\text{interact-abs}(Mat_i, Mat_j)$ is defined, from the conservative property of the abstract interact-abs function, we have $\text{interact}(\sigma(x_1), \sigma(x_2), r) \in \text{interact-abs}(Mat_i, Mat_j)$. Thus, from [11], we have $\text{interact}(\sigma(x_1), \sigma(x_2), r) \in \Gamma(x)$. From this and the consistency assumption for σ , we have $\text{consistent}(\Gamma, X \cup \{x\}, \sigma[x \mapsto \text{interact}(\sigma(x_1), \sigma(x_2), r)])$. Therefore, we have $\text{consistent}(\Gamma, X \setminus \{x_1, x_2\} \cup \{x\}, \sigma \setminus \{x_1, x_2\} [x \mapsto \text{interact}(\sigma(x_1), \sigma(x_2), r)])$.

Case for the rule T-SPLIT:

(7) $i = (\langle x_1, \dots, x_n \rangle = \text{split } x \text{ into } n)$

(8) $\Gamma, X \vdash x : \cup \overline{Mat_i}$

(9) $\Gamma(x) \subseteq \Gamma(x_1), \dots, \Gamma(x) \subseteq \Gamma(x_n)$

(10) $X' = X \setminus \{x\} \cup \{x_1, \dots, x_n\}$

Case analysis on [2]: There is only one case. Case of the rule E-SPLIT: We already have the typing judgement for the remaining statement s'' in [5]. By Lemma 6 on [8], we have $\Gamma(x) = \cup \overline{Mat_i}$ and $x \in X$. From the consistency condition, we have $\sigma(x) \in \Gamma(x)$. From the type of the split function, we have $\text{split}(\sigma(x), n) \in \Gamma(x)$. From this and [9], we have that for every $j \in \{1..n\}$, $\text{split}(\sigma(x), n) \in \Gamma(x_j)$. From this and the consistency assumption for σ , we have $\text{consistent}(\Gamma, X \cup \{x_1, \dots, x_n\}, \sigma[x_i \mapsto \text{split}(\sigma(x), n)])$. Therefore, we have $\text{consistent}(\Gamma, X \setminus \{x\} \cup \{x_1, \dots, x_n\}, (\sigma \setminus \{x\})[x_i \mapsto \text{split}(\sigma(x), n)])$.

Case for the rule T-If

(7) $i = \text{if } t \text{ then } s_1 \text{ else } s_2$

(8) $\Gamma, X \vdash t : \mathbb{N}$

(9) $\Gamma, X \vdash s_1, X_1$

(10) $\Gamma, X \vdash s_2, X_2$

(11) $X' = X_1 \cap X_2$

Case analysis on [2]: There are three cases.

Case of the rule E-IfR:

By Lemma 12, the type of t is preserved for t' . Thus, the rule T-If is applied to the new if instruction. The assumed consistency condition is preserved since the store σ stays unchanged in the step.

Case of the rule E-IfTRUE:

From [11], we have $X_1 \subseteq X'$. From this and [9], by Lemma 1, we have $\Gamma, X \vdash s_1, X'$. From this and [5], by Lemma 3, we have $\Gamma, X \vdash s_1 \bullet s'', X''$. The assumed consistency condition is preserved since the store σ stays unchanged in the step.

Case of the rule E-IfFALSE:

Similar to the previous case.

Case for the rule T-WHILE:

(7) $i = \text{while } t \text{ } s'''$

(8) $\Gamma, X \vdash \text{while } t \text{ } s''', X$

(9) $\Gamma, X \vdash t : \mathbb{N}$

(10) $\Gamma, X \vdash s''', X'''$

(11) $X \subseteq X'''$

(12) $X' = X$

Case analysis on [2]: There is only one case.

Case for the rule E-WHILE:

(13) $s' = \text{if } t \text{ then } (s''' \bullet \text{while } t \text{ } s'''; s'') \text{ else } s''$.

By Lemma 4 on [10], [11] and [8], we have

(14) $\Gamma, X \vdash (s''' \bullet \text{while } t \text{ } s'''), X_1$

(15) $X \subseteq X_1$

From [5] and [12], we have

(16) $\Gamma, X \vdash s'', X''$

By Lemma 2 on [14], [15] and [16], we have

(17) $\Gamma, X \vdash (s''' \bullet \text{while } t \text{ } s'''; s''), X_2$

(18) $X'' \subseteq X_2$

By the rule T-If on [9], [17], [16], we have

(19) $\Gamma, X \vdash \text{if } t \text{ then } (s''' \bullet \text{while } t \text{ } s'''; s'') \text{ else } s'', X_2 \cap X''$

From [19], [13] and [18], we have

$\Gamma, X \vdash s', X''$

The assumed consistency condition [3] is preserved since the store σ stays unchanged in the step.

□

2.4 Proof of Soundness

LEMMA 14 (SOUNDNESS OF TYPE INFERENCE FOR TERMS).

For every Γ, X, t, T, C and m ,
if

$\Gamma, X \vdash t : T \mid C$

m is a model for C

then

$m(\Gamma), X \vdash t : m(T)$

PROOF.

Hypothesis:

(1) $\Gamma, X \vdash t : T \mid C$

(2) m is a model for C

Structural induction on [1]:

Case the rule CT-VAR:

Trivial

Case the rule CT-MATH:

(3) $\Gamma, X \vdash t_1 : T_1 \mid C_1$

(4) $\Gamma, X \vdash t_2 : T_2 \mid C_2$

(5) $C = C_1 \cup C_2 \cup \{T_1 = T_2 = \mathbb{N} \vee T_1 = T_2 = \mathbb{R}\}$

(6) $T = T_1$

(7) $t = t_1 \oplus t_2$

From [2] and [5]:

(8) m is a model for C_1 .

(9) m is a model for C_2 .

By I.H. on ([3], [8]), ([4], [9]):

(10) $m(\Gamma), X \vdash t_1 : m(T_1)$

(11) $m(\Gamma), X \vdash t_2 : m(T_2)$

From [5], [2]:

$m(T_1) = m(T_2) = \mathbb{N} \vee m(T_1) = m(T_2) = \mathbb{R}$

We consider the first disjunct. The case for the second one is similar.

(12) $m(T_1) = m(T_2) = \mathbb{N}$

From [12], [10] and [11]:

(13) $m(\Gamma), X \vdash t_1 : \mathbb{N} \wedge m(\Gamma), X \vdash t_2 : \mathbb{N}$

By the rule T-MATH on [13], [14]:

(15) $m(\Gamma), X \vdash t_1 \oplus t_2 : \mathbb{N}$

From [15], [6], [7] and [12]:

(16) $m(\Gamma), X \vdash t : m(T)$

Case the rule CT-DETECT:

Similar to the previous case. The only interesting step is that considering the syntax of types T , the equality $m(T_1) \cap \{\mathbb{R}, \mathbb{N}\} = \emptyset$ implies $m(T_1) = \cup \overline{Mat}_i$.

1177 Case the rule CT-MAT:
1178 Trivial

1180 Case the rule CT-REAL:
1181 Trivial

1183 Case the rule CT-NAT:
1184 Trivial

1186 □

1187 **LEMMA 15 (SOUNDNESS OF TYPE INFERENCE FOR STATEMENTS).**

1188 *For every Γ, X, s, C, X'' , and m ,*

1189 *if*

1190 $\Gamma, X \vdash s, X'' \mid C$

1191 m is a model for C

1192 *then*

1193 $m(\Gamma), X \vdash s, X''$

1195 **PROOF.**

1196 Hypothesis:

1197 (1) $\Gamma, X \vdash s \mid C$

1198 (2) m is a model for C

1200 Structural induction on [1]:

1201 Case the rule CT-SKIP:

1202 Trivial.

1204 Case the rule CT-INST:

1205 (3) $s = i; s'$

1206 (4) $\Gamma, X \vdash i, X' \mid C_1$

1207 (5) $\Gamma, X' \vdash s', X'' \mid C_2$

1208 (6) $C = C_1 \cup C_2$

1209 From [2] and [6]:

1210 (7) m is a model for C_1

1211 (8) m is a model for C_2

1212 By I.H. on [5]:

1213 (9) $m(\Gamma), X' \vdash s, X''$

1214 We need to show

1215 (10) $m(\Gamma), X \vdash i, X'$

1216 and then by the rule T-INST on [9], [10]:

1217 (11) $m(\Gamma), X \vdash i; s', X''$

1218 and then by from [10], [11]:

1219 $m(\Gamma), X \vdash s, X''$

1221 Now assuming [4] and [7] that is

1222 (4) $\Gamma, X \vdash i, X' \mid C_1$

1223 (7) m is a model for C_1

1224 we show that

1225

1226 $m(\Gamma), X \vdash i, X'$

1227
1228 Case analysis on [4]:

1229
1230 Case the rule CT-ASSIGN-1:

1231 (12) $i = (x := v)$

1232 (13) $x : T \in \Gamma$

1233 (14) $\Gamma, X \vdash v : T' \mid C'$

1234 (15) $X' = X \cup \{x\}$

1235 (16) $C_1 = C' \cup \{T' \subseteq T\}$

1236 From [7] and [16]:

1237 (17) m is a model for C' .

1238 (18) $m(T') \subseteq m(T)$

1239 From [13]:

1240 (19) $x : m(T) \in m(\Gamma)$

1241 From Lemma 14 on [14] and [17]:

1242 (20) $m(\Gamma), X \vdash mat : m(T')$

1243 By the rule T-ASSIGN on [19], [20], and [18]:

1244 (21) $m(\Gamma), X \vdash x := mat, X \cup \{x\}$

1245 From [21], [15] and [12]:

1246 $m(\Gamma), X \vdash i, X'$

1247
1248 Case the rule CT-ASSIGN-2:

1249
1250 Similar to the previous case. The the rule T-ASSIGN-2 is used.

1251
1252 Case the rule CT-ASSIGN-3:

1253
1254 Similar to the case the rule CT-ASSIGN-1. To use the rule T-ASSIGN, in contrast to the case
1255 the rule CT-ASSIGN-1 that proved $t = v$, the disjunct $T' = \mathbb{R} \vee T' = \mathbb{N}$ is proved to use the
1256 rule T-ASSIGN-3.

1257
1258 Case the rule CT-Mix:

1259 (12) $\Gamma, X \vdash x_1 : T \mid C$

1260 (13) $\Gamma, X \vdash x_2 : T' \mid C'$

1261 (14) $\Gamma, X \vdash t : T'' \mid C''$

1262 (15) $i = (x := \text{mix } x_1 \text{ with } x_2 \text{ for } t)$

1263 (16) $X' = X \setminus \{x_1, x_2\} \cup \{x\}$

1264 (17) $C_1 = C \cup C' \cup C'' \cup \{T \cap \{\mathbb{R}, \mathbb{N}\} = \emptyset, T' \cap \{\mathbb{R}, \mathbb{N}\} = \emptyset, T'' = \mathbb{R}$

1265 $\frac{Mat_i \in T \wedge Mat_j \in T' \Rightarrow \text{interact-abs}(Mat_i, Mat_j) \subseteq \Gamma(x)}{}$

1266 From [7] and [17]:

1267 (18) m is a model for C .

1268 (19) m is a model for C' .

1269 (20) m is a model for C'' .

1270 (21) $m(T) \cap \{\mathbb{R}, \mathbb{N}\} = \emptyset$

1271 (22) $m(T') \cap \{\mathbb{R}, \mathbb{N}\} = \emptyset$

1272 (23) $m(T'') = \mathbb{R}$

$$(24) \overline{Mat_i \in m(T) \wedge Mat_j \in m(T') \Rightarrow \text{interact-abs}(Mat_i, Mat_j) \subseteq m(\Gamma(x))}$$

From Lemma 14 on [12] and [18]:

$$(25) m(\Gamma), X \vdash x_1 : m(T)$$

From Lemma 14 on [13] and [19]:

$$(26) m(\Gamma), X \vdash x_2 : m(T')$$

From Lemma 14 on [14] and [20]:

$$(27) m(\Gamma), X \vdash t : m(T'')$$

From [21], there exists i such that:

$$(28) m(T) = \cup \overline{Mat_i}$$

From [22], there exists j such that:

$$(29) m(T') = \cup \overline{Mat_j}$$

From [25] and [28]:

$$(30) m(\Gamma), X \vdash x_1 : \cup \overline{Mat_i}$$

From [26] and [29]:

$$(31) m(\Gamma), X \vdash x_2 : \cup \overline{Mat_j}$$

From [27] and [23]:

$$(32) m(\Gamma), X \vdash t : \mathbb{R}$$

From [24], [28], and [29]:

$$(33) \text{interact-abs}(Mat_i, Mat_j) \subseteq \Gamma(x)$$

By the rule T-Mix on [30], [31], [32], and [33]:

$$(34) \Gamma, X \vdash x := \text{mix } x_1 \text{ with } x_2 \text{ for } t, X \setminus \{x_1, x_2\} \cup \{x\}$$

From [34], [15] and [16]:

$$(34) \Gamma, X \vdash i, X'$$

Case the rule CT-SPLIT:

Similar to the case for the rule CT-Mix. The Lemma 14 is applied to the constraint typing judgement for x .

Case the rule CT-IF:

Immediate from applying the Lemma 14 on t and the induction hypothesis on s_1 and s_2 .

Case the rule CT-WHILE:

Immediate from applying the Lemma 14 on t and the induction hypothesis on s .

□

2.5 Proof of Completeness

LEMMA 16 (COMPLETENESS OF TYPE INFERENCE FOR TERMS).

For all Γ, X, t, T , and C ,

if

$$m(\Gamma), X \vdash t : T$$

$$\Gamma, X \vdash t : T' \mid C$$

then

m is a model for C

$$m(T') = T$$

PROOF.

Hypothesis

$$(1) m(\Gamma), X \vdash t : T$$

$$(2) \Gamma, X \vdash t : T' \mid C$$

Proof by induction on the given constraint typing derivation [2].

Case for the rule CT-VAR:

$$(3) t = x$$

$$(4) x : T \in \Gamma$$

$$(5) x \in X$$

$$(6) C = \emptyset$$

By the inversion Lemma 6 on [1]:

$$(7) x : T' \in m(\Gamma)$$

$$(8) x \in X$$

From [4] and [7]:

$$(7) T' = m(T)$$

The conclusion is immediate from [6] and [7].

Case for the rule CT-MATH:

$$(3) \Gamma, X \vdash t_1 : T_1 \mid C_1$$

$$(4) \Gamma, X \vdash t_2 : T_2 \mid C_2$$

$$(5) C = C_1 \cup C_2 \cup \{T_1 = T_2 = \mathbb{N} \vee T_1 = T_2 = \mathbb{R}\}$$

$$(6) T' = T_1$$

By the inversion Lemma 6 on [1]:

$$(7) m(\Gamma), X \vdash t_1 : T$$

$$(8) m(\Gamma), X \vdash t_2 : T$$

$$T = \mathbb{N} \vee T = \mathbb{R}$$

We consider the case for the first disjunct:

(The case for the second one is similar.)

$$(9) T = \mathbb{N}$$

By induction hypothesis on [7] and [3]:

$$(10) m \text{ is a model for } C_1$$

$$(11) m(T_1) = T$$

By induction hypothesis on on [8] and [4]:

$$(12) m \text{ is a model for } C_2$$

$$(13) m(T_2) = T$$

From [11], [13] and [9]:

$$(14) m(T_1) = m(T_2) = \mathbb{N}$$

From [5], [10], [12] and [14]:

$$(15) m \text{ is a model for } C.$$

From [6], [11], and [9]:

$$(16) m(T') = T$$

The conclusion is [15] and [16].

Case for the rule CT-DETECT:

Similar to the case for the rule T-MATH. By induction hypothesis on x and t .

Case for the rule CT-MAT

Trivial.

Case for the rule CT-REAL:

Trivial.

Case for the rule CT-NAT:

Trivial.

□

LEMMA 17 (COMPLETENESS OF TYPE INFERENCE FOR STATEMENTS).

For all Γ, X, s, X', C , and m ,

if

$m(\Gamma), X \vdash s, X'$

$\Gamma, X \vdash s, X'' \mid C$

then

m is a model for C .

PROOF.

Hypothesis

(1) $m(\Gamma), X \vdash s, X_1$

(2) $\Gamma, X \vdash s, X_2 \mid C$

We show that

m is a model for C .

Induction on the derivation of [2]:

Case for the rule CT-SKIP:

Trivial.

Case for the rule CT-INST:

(3) $s = i; s'$

(4) $\Gamma, X \vdash i, X' \mid C_1$

(5) $\Gamma, X' \vdash s', X_2 \mid C_2$

(6) $C = C_1 \cup C_2$

From [1] and [3]:

(7) $m(\Gamma), X \vdash i; s', X_1$

By the inversion Lemma 7 on [7]:

(8) $m(\Gamma), X \vdash i, X''$

(9) $m(\Gamma), X'' \vdash s', X_2$

By induction hypothesis on [9] and [5]:

(10) m is a model for C_2

We will show that from [8] and [4]:

(11) m is a model for C_1

From [6], [10] and [11]:

m is a model for C

Hypothesis:

(8) $m(\Gamma), X \vdash i, X''$

(4) $\Gamma, X \vdash i, X' \mid C_1$

Conclusion:

m is a model for C_1

Case analysis on derivation of [4]:

Case for the rule CT-Assign-1:

(12) $i = (x := v)$

(13) $x : T \in \Gamma$

(14) $\Gamma, X \vdash v : T' \mid C'$

(15) $X' = X \cup \{x\}$

(16) $C_1 = C' \cup \{T' \subseteq T\}$

From [8] and [12]:

(17) $m(\Gamma), X \vdash x := v, X''$

By the inversion Lemma 7 on [17]:

(18) $x : T'' \in m(\Gamma)$

(19) $m(\Gamma), X \vdash v : T'''$

(22) $T''' \subseteq T''$

(23) $X'' = X \cup \{x\}$

From [18], [13]:

(24) $m(T) = T''$

By Lemma 16 on [19], [14]:

(25) m is a model for C'

(26) $m(T') = T'''$

From [22], [24] and [26]

(27) $m(T') \subseteq m(T)$

From [16], [24], [25] and [27]

m is a model for C_1

Case for the rule CT-Assign-2:

Similar to the case for the rule CT-Assign-1.

Case for the rule CT-Assign-3:

Similar to the case for the rule CT-Assign-1.

Case for the rule CT-Mix:

(12) $i = (x := \text{mix } x_1 \text{ with } x_2 \text{ for } t)$

(12) $\Gamma, X \vdash x_1 : T \mid C$

(13) $\Gamma, X \vdash x_2 : T' \mid C'$

(14) $\Gamma, X \vdash t : T'' \mid C''$

$$(15) \frac{C_1 = C \cup C' \cup C'' \cup \{T \cap \{\mathbb{R}, \mathbb{N}\} = \emptyset, T' \cap \{\mathbb{R}, \mathbb{N}\} = \emptyset, T'' = \mathbb{R}\} \quad \text{Mat}_k \in T \wedge \text{Mat}_l \in T' \Rightarrow \text{interact-abs}(\text{Mat}_k, \text{Mat}_l) \subseteq \Gamma(x)}{}$$

From [8] and [12]:

$$(16) m(\Gamma), X \vdash x := \text{mix } x_1 \text{ with } x_2 \text{ for } t, X''$$

By the inversion Lemma 7 on [16]:

$$(17) m(\Gamma), X \vdash x_1 : \overline{\cup \text{Mat}_i}$$

$$(18) m(\Gamma), X \vdash x_2 : \overline{\cup \text{Mat}_j}$$

$$(19) m(\Gamma), X \vdash t : \mathbb{R}$$

$$(20) \overline{\text{interact-abs}(\text{Mat}_i, \text{Mat}_j) \subseteq m(\Gamma(x))}$$

By Lemma 16 on [17], [12]:

$$(21) m \text{ is a model for } C$$

$$(22) m(T) = \overline{\cup \text{Mat}_i}$$

By Lemma 16 on [18], [13]:

$$(23) m \text{ is a model for } C'$$

$$(24) m(T') = \overline{\cup \text{Mat}_j}$$

By Lemma 16 on [19], [14]:

$$(25) m \text{ is a model for } C''$$

$$(26) m(T'') = \mathbb{R}$$

From [22], [24], and [20]:

$$(27) \text{Mat}_k \in m(T) \wedge \text{Mat}_l \in m(T') \Rightarrow \text{interact-abs}(\text{Mat}_k, \text{Mat}_l) \subseteq m(\Gamma(x))\}$$

From [22], [24], and [26]:

$$(28) m(T) \cap \{\mathbb{R}, \mathbb{N}\} = \emptyset, m(T') \cap \{\mathbb{R}, \mathbb{N}\} = \emptyset, m(T'') = \mathbb{R}$$

From [15], [21], [23], [25], [28] and [27]:

$$m \text{ is a model for } C_1$$

Case for the rule CT-SPLIT:

Similar to the case for the rule CT-MIX.

Case for the rule CT-IF:

Immediate from the induction hypothesis on s_1 and s_2 and Lemma 16 on t .

Case for the rule CT-WHILE:

Immediate from the induction hypothesis on s and Lemma 16 on t .

□

3 SYNTAX STUDY: ELISA PROTOCOLS

The *enzyme-linked immunosorbent assay (ELISA)* is a test that uses antibodies and color changes to identify a substance. ELISA assays are commonly performed by interacting a group of chemicals with an immobile antibody to detect the presence and concentration of particular opiates. A technique was introduced, to allow these assays to be performed on DMFB technology by baking the enzyme directly onto the top plate of the DMFB, requiring further syntactic extensions to declare stationary substances and to explicitly move mobile substances to interact with them.

8 shows an example ELISA assay, which is one of the steps of the larger hierarchical opiate-detection decision tree. The operator will need to swap out top plates with different antigens for each ELISA assay, or use a very large electrowetting array to support a top plate to which all necessary antigens for all of the ELISA assays are affixed.

Biologists use many equivalent terms that mean the same thing as Mix. BioScript supports several of these terms, including Tap and Vortex, as shown in 8. The execution engine recognizes these terms as being equivalent and converts them all to its own internal mix operation.

The repeat operation is part of BioScript's core instruction set. It allows the programmer to specify a sequence of instructions that will be repeated a constant number of times, which is far more common in biochemistry than in general computer programming. Thus, instead of the complex, unconstrained loop syntax employed by modern programming languages, we have opted for a simple syntax to make BioScript accessible to scientists with limited programming experience.

The Incubate and Heat operations leverage the external heaters that are integrated into or placed near the DMFB. A DMFB that lacks heating capabilities cannot perform this ELISA assay. Thus, the compiler writer must specialize a language syntax for each DMFB variant based on its integrated peripherals.

The final operation in any ELISA is the detection phase that compares the reading obtained from the sample to a reading obtained from a control. This provides a certain measurement of the presence and concentration of the opiate in the sample. In the hierarchical decision tree immunoassay, the result of this comparison determines which ELISA assay to execute next.

Stationary: Anti-Fentanyl

Move 20uL of Urine Sample to Anti-Fentanyl

Move 100uL of Fentanyl-Conjugate to Anti-Fentanyl

Tap Anti-Fentanyl for 60s

Incubate Anti-Fentanyl at 23°C for 60min

Drain Anti-Fentanyl

Repeat 6 times {

Move 350uL of Distilled Water to Anti-Fentanyl

Vortex Anti-Fentanyl for 45s

Drain Anti-Fentanyl

}

Move 100uL of TMB Substrate to Anti-Fentanyl

Incubate Anti-Fentanyl at 23C for 30min

Mix Anti-Fentanyl with 100uL of Stop Reagent

for 60s

Negative Reading = Measure the fluorescence of

Anti-Fentanyl for 30min

Drain Anti-Fentanyl

Fig. 8. Fentanyl ELISA




Fig. 9. Depiction of the electrowetting principle [29, 37]: applying an electrostatic potential to a droplet at rest reduces the contact angle with the surface, thereby increasing the surface area in contact with the droplet.




Fig. 10. Left: A digital microfluidic biochip (DMFB) is a planar array of electrodes [15, 19, 36, 39, 41]. Right: Cross-sectional view.




Fig. 11. A droplet is transported from CE2 to CE3 by activating CE3, and then deactivating CE2 (white: activated electrode; black: deactivated electrode).




Fig. 12. Five basic fluidic operations.

4 BACKGROUND

4.1 Digital Microfluidic Technology

As a proof of concept, the current version of the BioScript compiler targets *digital microfluidic biochips (DMFBs)*, which are LoCs based on *electrowetting on dielectric (EWoD)* technology. Fig. 9 illustrates the electrowetting effect: applying an electrical potential to a liquid droplet resting on a hydrophobic surface decreases the contact angle with the surface. This increases the "wet" area of contact between the droplet and the surface [29, 37]. As shown in Fig. 10, a DMFB is a planar array of electrodes which leverages the electrowetting principle to induce droplet transport [15, 19, 36, 39, 41]. Activating and deactivating nearby electrodes in sequence induces droplet transport, as shown in Fig. 11.

A DMFB has a small *instruction set architecture (ISA)* consisting of five primitive operations, as shown in Fig. 12: droplet transport, splitting and merging two droplets, actively mixing two (or more) droplets by first merging them and moving them (typically in a confined space around a pivot), and droplet storage. Droplet input from reservoirs on the perimeter of the chip and output (including waste/disposal) can also be controlled electrically. Since the operations in Fig. 12 can be performed anywhere on the DMFB, a DMFB can be interpreted as a reconfigurable spatial computing device, similar in principle to an FPGA, as the function performed by a (group) of electrode(s) may change over time.

Additional operations can be realized by integrating external devices such as heaters [32], optical detectors [31], and a variety of integrated sensors [27, 38, 45, 46]. In particular, the hierarchical




Fig. 13. Illustration of a typical DMFB compiler for assays without control flow (DAGs). The compiler must schedule, place, and route all operations to produce an executable sequence of electrode activations to automatically execute the assay (activated electrodes typically have a droplet on top).




Fig. 14. A dynamic execution model for cyber-physical DMFBs.

decision tree opiate detection immunoassay will require optical sensors for *immunofluorescence* and impedimetric immunosensors utilizing *electrochemical impedance spectroscopy (EIS)* [5, 10, 51].

Despite the lack of high-level language support, compilers for assays without control flow are mature. As shown in Fig. 13 the compiler must solve three interdependent NP-complete problems: operation scheduling [11, 16, 30, 40, 42, 49], reconfigurable module placement [8, 17, 28, 33, 34, 48, 52, 53, 55], and droplet routing [7, 9, 22, 25, 26, 43, 44, 50, 56]. The scheduler must determine the time steps at which each biochemical operation occurs, while satisfying droplet dependency constraints and physical resource constraints of the device. The placer determines the location on the 2D electrode array where each operation is performed as the reaction progresses over time. The router ensures that droplets are transported from their start/stop points at appropriate times during execution of the chemical reaction, while ensuring that droplets undergoing transport do not inadvertently collide with one another or any other ongoing operations on the chip. The router may also introduce wash droplets to remove residue left by droplets that travel over the surface of the chip, to prevent cross-contamination [23, 54, 57].

A *cyber-physical* DMFB leverages integrated sensors [31, 38, 45, 46] and/or real-time video monitoring [6, 13, 21, 47] to send online feedback in real-time to a PC controller. Cyber-physical integration provides the capability of control flow [18]: the assay designer can specify different steps to take based on the sensor feedback.

As control flow cannot be predicted statically, the simple compiler model shown in Fig. 13 no longer suffices. In response, a software-based dynamic execution model, shown in Fig. 14 was introduced [18]. The assay, which is now specified as a *control flow graph* is compiled and executed one basic block at a time. Each basic block other than the CFG exit node terminates with a condition, typically based on sensory feedback. The control software evaluates the condition, determines the next basic block to execute, quickly compiles it on-the-fly, and then initiates execution.

5 ASSAY EXECUTION VIDEOS

Below is a list of the videos generated while gathering data. We are not claiming these as novel or interesting contributions. They are merely provided to demonstrate how the simulator executes a given assay.

Urine Opiate Panel with initial positive indication: https://1drv.ms/v/s!AqwZF7jQGx_IgahNYojptiZr4RpjOg
Urine Opiate Panel with initial negative indication: https://1drv.ms/v/s!AqwZF7jQGx_IhJRtiOKGuSTm0qXX0
PCR Assay: https://1drv.ms/v/s!AqwZF7jQGx_IhNgWWxzt8GwtfICcFw
Image Probe Synthesis: https://1drv.ms/v/s!AqwZF7jQGx_IhNgbKsUIzDGqfyOhbw
Neurotransmitter Sensing: https://1drv.ms/v/s!AqwZF7jQGx_IhNgZjBF4vlsnDr1uCA
PCR Droplet Replenishment Assay: https://1drv.ms/v/s!AqwZF7jQGx_IhJRqfq1dLhbBUIL9RQ
Probabilistic PCR Full: https://1drv.ms/v/s!AqwZF7jQGx_IhJRo-xDbVT_Q6UsXPg
Probabilistic PCR Early Exit: https://1drv.ms/v/s!AqwZF7jQGx_IhJRpLXveGI8Qr7BWrw
Broad Spectrum Opiate Panel: https://1drv.ms/v/s!AqwZF7jQGx_IhJR64f6laDx16T_QJQ
Fentanyl ELISA Assay: https://1drv.ms/v/s!AqwZF7jQGx_IgahJFK519Id1HCad4g
Ciprofloxacin ELISA Assay: https://1drv.ms/v/s!AqwZF7jQGx_IgahKrD6EDWPbsy9G0g
Heroin ELISA Assay: https://1drv.ms/v/s!AqwZF7jQGx_IgahIrCW0z08KR1111A
Morphine ELISA Assay: https://1drv.ms/v/s!AqwZF7jQGx_IgahG3PHVTctW7u0rOw
Oxycodone ELISA Assay: https://1drv.ms/v/s!AqwZF7jQGx_IgahHIgmGgR43A-17EA

6 BIOSCRIPT'S CHEMTYPE TESTS

6.1 Real-world Prevention Tests

AIHA Example1 where ChemType would have prevented an incident

Tetrachloroethylene = 50uL of C(=C(Cl)Cl)(Cl)Cl

Nitric_Acid = 50uL of [N+](=O)(O)[O-]

a = Mix Tetrachloroethylene with Nitric_Acid for 40s

AIHA Example2 where ChemType would have prevented an incident

Methanol = 50uL of CO

Nitric_Acid = 50uL of [N+](=O)(O)[O-]

a = Mix Methanol with Nitric_Acid for 40s

AIHA Example3 where ChemType would have prevented an incident

Nitrogen = N#N

Diaminopropane = 50uL of CC(CN)N

Potassium_Hydride = 50uL of [H-].[K+]

a = Mix Diaminopropane and Potassium_Hydride for 20s

b = Mix a and Nitrogen for 10s

Real world example where ChemType would have prevented an incident

Dichlor = 50uL CC(=CCl)Cl

Calcium_Hypo = 50uL of [O-]Cl.[O-]Cl.[Ca+2]

```

1716
1717 Mix Dichlor and Calcium_Hypo for 20s
1718 =====
1719
1720 6.2 Synthetic Failures
1721 Example where ChemType will not allow the mix to occur
1722 Nitric_Acid = 50uL of [N+](=O)(O)[O-]
1723 Hydrochloric_Acid = 50uL of Cl
1724
1725 a = Mix Nitric_Acid with Hydrochloric_Acid for 40s
1726 =====
1727 Example where ChemType will not allow the mix to occur
1728 Isopropanol = 50uL of CC(C)O
1729 Hydrochloric_Acid = 50uL of Cl
1730
1731 a = Mix Isopropanol with Hydrochloric_Acid for 40s
1732 =====
1733 Example where ChemType will not allow the mix to occur
1734 Formaldehyde = 50uL of C=O
1735 Formaldehyde2 = 50uL of C=O
1736
1737 a = Mix Formaldehyde with Formaldehyde2 for 40s
1738 =====
1739 Example where ChemType will not allow the mix to occur
1740 Benzoin = 50uL of C1=CC=C(C=C1)C(C(=O)C2=CC=CC=C2)O
1741 Urea = 50uL of C(=O)(N)N
1742 Formaldehyde = 50uL of C=O
1743
1744 a = Mix Benzoin with Urea for 40s
1745 b = Mix Formaldehyde with a for 20s
1746 =====
1747 Example where ChemType will not allow the mix to occur
1748 Benzene = 50uL of C1=CC=CC=C1
1749 Formaldehyde = 50uL of C=O
1750 Formaldehyde2 = 50uL of C=O
1751
1752 a = Mix Benzene with Formaldehyde for 40s
1753 b = Mix Formaldehyde2 with a for 20s
1754 =====
1755 Example where ChemType will not allow the mix to occur
1756 Benzene = 50uL of C1=CC=CC=C1
1757 Urea = 50uL of C(=O)(N)N
1758 Hydrochloric Acid = 50uL of Cl
1759
1760 a = Mix Benzene with Urea for 40s
1761
1762
1763
1764

```

```
1765 b = Mix Hydrochloric Acid with a for 20s
1766 =====
1767     Example where ChemType would have prevented an incident
1768 Hydrogen_Peroxide = 50uL O2
1769 Sulfuric_Acid = 50uL of OS(=O)(=O)O
1770 Acetone = 50uL of CC(=O)C
1771
1772 a = Mix Hydrogen_Peroxide and Sulfuric_Acid for 20s
1773 b = Mix a and Acetone for 10s
1774 =====
1775
```

6.3 Synthetic Passes

```
1777 Example where ChemType will allow the mix to occur
1778 Benzene = 50uL of C1=CC=CC=C1
1779 Urea = 50uL of C(=O)(N)N
1780
1781 a = Mix Benzene with Urea for 40s
1782 =====
1783     Example where ChemType will allow the mix to occur
1784 Benzene = 50uL of C1=CC=CC=C1
1785 Urea = 50uL of C(=O)(N)N
1786 Benzotrichloride = 50uL of C1=CC=C(C=C1)C(Cl)(Cl)Cl
1787
1788 a = Mix Benzene with Urea for 40s
1789 b = Mix Benzotrichloride with a for 20s
1790 =====
1791
```

7 BIOSCRIPT ASSAYS

7.1 Syntactic Sugar Translation

For the purposes of adhering to scientific assay nomenclature, we have chosen to represent the following assays in the way a biological scientist might. Adhering to the community's "syntax" requires some syntactic sugar. This mapping is described in Table 1.

```
1799 Thus, given:
1800 move x to y
1801 tap y
1802 BioScript would condense into:
1803 y := mix x with y
1804
```

7.2 Urine Opiate Hierarchy

Two main paths based on outcome of Broad Spectrum Opiate(BSO) Test.

When(BSO) Test Indicates Positive: https://1drv.ms/v/s!AqwZF7jQGx_IgahNYojptiZr4RpjOg

When(BSO) Test Indicates Negative: https://1drv.ms/v/s!AqwZF7jQGx_IhJRtiOKGuSTm0qXX0g

General Instruction	BioScript Equivalent
stationary Anti-Morphine	Anti-Morphine1 := x µl of Anti-Morphine
heat x at n	Not yet modeling temperature\volume. Left for future work.
incubate x at n	Not yet modeling temperature\volume. Left for future work.
move x to y	y := mix x with y
tap x	x := mix x with x
invert x	x := mix x with x
vortex x	x := mix x with x
drain x	freshVariable := x
repeat x times	while x > 0
measure x of x'	detect x on x'
perform...	x := detect x on x'

Table 1. Syntactic sugar translations from assay to BioScript.

```
BioScript Code:
Stationary Anti-Morphine
Stationary Anti-Oxy
Stationary Anti-Fentanyl
Stationary Anti-Ciprofloxacin
Stationary Anti-Heroin
Stationary Oxycodone-Enzyme

Stop-Reagent1 = 10uL of Stop-Reagent
Stop-Reagent2 = 10uL of Stop-Reagent
Stop-Reagent3 = 10uL of Stop-Reagent
Stop-Reagent4 = 10uL of Stop-Reagent
Stop-Reagent5 = 10uL of Stop-Reagent

Move 10uL of UrineSample to Anti-Morphine
Move 10uL of UrineSample to Anti-Oxy
Move 10uL of UrineSample to Anti-Fentanyl
Move 10uL of UrineSample to Anti-Ciprofloxacin
Move 10uL of UrineSample to Anti-Heroin

MorphineReading = Measure the fluorescence of Anti-Morphine for 5s
OxyReading = Measure the fluorescence of Anti-Oxy for 5s
FentanylReading = Measure the fluorescence of Anti-Fentanyl for 5s
CiproReading = Measure the fluorescence of Anti-Ciprofloxacin for 5s
HeroinReading = Measure the fluorescence of Anti-Heroin for 5s

Drain Anti-Morphine
Drain Anti-Oxy
Drain Anti-Fentanyl
Drain Anti-Ciprofloxacin
Drain Anti-Heroin
```

```
1863
1864 if (((FentanylReading == true) or (MorphineReading == true)) or (((OxyReading == true) or (
1865
1866     Move 20uL of Urine-Sample to Anti-Heroin
1867     Move 100uL of Heroin-Conjugate to Anti-Heroin
1868     Tap Anti-Heroin for 60s
1869     Incubate Anti-Heroin at 23C for 60min
1870     Drain Anti-Heroin
1871
1872     Repeat 6 times {
1873         Move 350uL of Distilled-Water to Anti-Heroin
1874         Vortex Anti-Heroin for 45s
1875         Drain Anti-Heroin
1876     }
1877
1878     Move 100uL of TMB-Substrate to Anti-Heroin
1879     Incubate Anti-Heroin at 23C for 30min
1880     Mix Anti-Heroin with 100uL of Stop Reagent1 for 60s
1881
1882     Heroin Reading = Measure the fluorescence of Anti-Heroin for 30min
1883     Drain Anti-Heroin
1884
1885     Move 20uL of Urine-Sample to Oxycodone-Enzyme
1886     Move 100uL of Oxycodone-Conjugate to Oxycodone-Enzyme
1887     Tap Oxycodone-Enzyme for 60s
1888     Incubate Oxycodone-Enzyme at 23C for 60min
1889     Drain Oxycodone-Enzyme
1890
1891     Repeat 6 times {
1892         Move 350uL of Distilled-Water to Oxycodone-Enzyme
1893         Vortex Oxycodone-Enzyme for 45s
1894         Drain Oxycodone-Enzyme
1895     }
1896
1897     Move 100uL of TMB-Substrate to Oxycodone-Enzyme
1898     Incubate Oxycodone-Enzyme at 23C for 30min
1899     Mix Oxycodone-Enzyme with 100uL of Stop-Reagent2 for 60s
1900
1901     Oxy Reading = Measure the fluorescence of Oxycodone-Enzyme for 30min
1902     Drain Oxycodone-Enzyme
1903
1904     if ((Heroin Reading == false) and (Oxy Reading == false)) {
1905         Move 20uL of UrineSample to Anti-Ciprofloxacin
1906         Move 100uL of Ciprofloxacin-Conjugate to Anti-Ciprofloxacin
1907         Tap Anti-Ciprofloxacin for 60s
1908         Incubate Anti-Ciprofloxacin at 23C for 60min
1909         Drain Anti-Ciprofloxacin
1910
1911
```

```

1912 Repeat 5 times {
1913     Move 250uL of Distilled-Water to Anti-Ciprofloxacin
1914     Vortex Anti-Ciprofloxacin for 45s
1915     Drain Anti-Ciprofloxacin
1916 }
1917
1918 Move 50uL of TMB-Substrate to Anti-Ciprofloxacin
1919 Incubate Anti-Ciprofloxacin at 25C for 30min
1920 Mix Anti-Ciprofloxacin with 100L of Stop Reagent3 for 60s
1921
1922 Urine Reading = Measure the fluorescence of Anti-Ciprofloxacin for 5min
1923 Drain Anti-Ciprofloxacin
1924 }
1925
1926 } else {
1927     Split Epithelial = Mix UrineSample with Azoxymethane for 5s
1928     Move Split Epithelial to Anti-Fentanyl
1929     Move 100uL of Fentanyl-Conjugate to Anti-Fentanyl
1930     Tap Anti-Fentanyl for 60s
1931     Incubate Anti-Fentanyl at 23C for 60min
1932     Drain Anti-Fentanyl
1933
1934     Repeat 6 times {
1935         Move 350uL of Distilled-Water to Anti-Fentanyl
1936         Vortex Anti-Fentanyl for 45s
1937         Drain Anti-Fentanyl
1938     }
1939
1940     Move 100uL of TMB-Substrate to Anti-Fentanyl
1941     Incubate Anti-Fentanyl at 23C for 30min
1942     Mix Anti-Fentanyl with 100uL of Stop Reagent4 for 60s
1943     Urine Reading = Measure the fluorescence of Anti-Fentanyl for 30min
1944     Drain Anti-Fentanyl
1945
1946     if (!(Heroin Reading >= FentanylControl)) {
1947         Move 20uL of Urine-Sample to Oxycodone-Enzyme
1948         Move 100uL of Oxycodone-Conjugate to Oxycodone-Enzyme
1949         Tap Oxycodone-Enzyme for 60s
1950         Incubate Oxycodone-Enzyme at 23C for 60min
1951         Drain Oxycodone-Enzyme
1952
1953         Repeat 6 times {
1954             Move 350uL of Distilled-Water to Oxycodone-Enzyme
1955             Vortex Oxycodone-Enzyme for 45s
1956             Drain Oxycodone-Enzyme
1957         }
1958
1959         Move 100uL of TMB-Substrate to Oxycodone-Enzyme
1960

```



```

1961     Incubate Oxycodone-Enzyme at 23C for 30min
1962     Mix Oxycodone-Enzyme with 100uL of Stop Reagent5 for 60s
1963
1964     Urine Reading = Measure the fluorescence of Oxycodone-Enzyme for 30min
1965     Drain Oxycodone-Enzyme
1966 }
1967 }
1968

```

1969 7.3 PCR Droplet Replenishment

1970 PCR Droplet Replenishment Assay: https://1drv.ms/v/s!AqwZF7jQGx_IhJRqfq1dLhbBUIL9RQ
 1971 BioScript Code:

```

1972 PCRMix = Vortex 50 uL of PCRMasterMix with 50 uL of Template for 1s
1973
1974 Repeat 50 times {
1975     Heat PCRMix at 95C for 20s
1976     volumeWeight = Weigh PCRMix
1977
1978     if (volumeWeight <= 50uL) {
1979         replacement = Vortex 25uL of PCRMasterMix with 25uL of Template for 5s
1980         Heat replacement at 95C for 45s
1981         PCRMix = Mix PCRMix with replacement for 5s
1982     }
1983
1984     Heat PCRMix at 68C for 30s
1985     Heat PCRMix at 95C for 45s
1986 }
1987
1988 Heat PCRMix at 68C for 5min
1989 Save PCRMix
1990

```

1991 7.4 Probabilistic PCR

```

1992 Probabilistic PCR Full: https://1drv.ms/v/s!AqwZF7jQGx\_IhJRo-xDbVT\_Q6UsXPg
1993 Probabilistic PCR Early Exit: https://1drv.ms/v/s!AqwZF7jQGx\_IhJRpLXveGI8Qr7BWrw
1994
1995 PCR-Master-Mix = Mix 50uL of PCRMix with 50uL of Buffer
1996
1997 # Initialization Step
1998 Heat PCR-Master-Mix at 94C for 2min
1999
2000 # Denaturation
2001 Repeat 20 times {
2002     Heat PCR-Master-Mix at 94C for 20s
2003     Heat PCR-Master-Mix at 50C for 40s
2004 }
2005 DNA Sensor = Measure the fluorescence of PCR-Master-Mix for 30s
2006 if (DNA Sensor <= 85) {
2007     Drain PCR-Master-Mix
2008 }
2009

```

```
2010
2011 Repeat 20 times {
2012     Heat PCR-Master-Mix at 94C for 20s
2013     Heat PCR-Master-Mix at 50C for 40s
2014 }
2015 # Final Elongation
2016 Heat PCR-Master-Mix at 70C for 5min
2017
2018 Save PCR-Master-Mix
2019
2020 7.5 ChemType Test1
2021 Allyl Ethyl Ether = 50uL of C=CCOC1=CC=CC=C1
2022
2023 Cycloheptatriene = 50uL of C1C=CC=CC=C1
2024
2025 Cycloheptatriene2 = 50uL of C1C=CC=CC=C1
2026
2027 Sulfuric Acid = 50uL of OS(=O)(=O)O
2028
2029 Chlorotrifluoromethane = 50uL of C(F)(F)(F)Cl
2030
2031 Ammonium dichromate = 50uL of [NH4+].[NH4+].[O-][Cr](=O)(=O)O[Cr](=O)(=O)[O-]
2032
2033
2034 a = Mix Allyl Ethyl Ether with Cycloheptatriene for 50s
2035 b = Vortex a with Sulfuric Acid for 30s
2036 c = Vortex b with Chlorotrifluoromethane for 45s
2037 d = Mix Cycloheptatriene2 with Ammonium dichromate
2038
2039 7.6 ChemType Test2
2040 Isopropyl-Alcohol = 50uL of CC(C)O
2041
2042 Isopropyl-Alcohol2 = 50uL of CC(C)O
2043
2044 Allyl-Ethyl-Ether = 50uL of C=CCOC1=CC=CC=C1
2045
2046 Cycloheptatriene = 50uL of C1C=CC=CC=C1
2047
2048 Cycloheptatriene2 = 50uL of C1C=CC=CC=C1
2049
2050 Cycloheptatriene3 = 50uL of C1C=CC=CC=C1
2051
2052 Sulfuric-Acid = 50uL of OS(=O)(=O)O
2053
2054 Chlorotrifluoromethane = 50uL of C(F)(F)(F)Cl
2055
2056 Ammonium-dichromate = 50uL of [NH4+].[NH4+].[O-][Cr](=O)(=O)O[Cr](=O)(=O)[O-]
2057
2058
```

2059
2060 a = Mix Cycloheptatriene with Isopropyl-Alcohol for 30s
2061 a2 = Mix Cycloheptatriene2 with Isopropyl-Alcohol2 for 30s
2062
2063 b = Vortex Cycloheptatriene3 with Chlorotrifluoromethane
2064 c = Tap a with b
2065 d = Invert c with Allyl-Ethyl-Ether
2066 e = Mix a2 with Sufluric-Acid
2067

2068 7.7 ChemType Test3

2069 Sodium-Hydroxide = 50uL of [OH-].[Na+]
2070
2071 Isopropyl-Alcohol = 50uL of CC(C)O
2072
2073 Allyl-Ethyl-Ether = 50uL of C=CCOC1=CC=CC=C1
2074
2075 Cycloheptatriene = 50uL of C1C=CC=CC=C1
2076
2077 Cycloheptatriene2 = 50uL of C1C=CC=CC=C1
2078
2079 Sufluric-Acid = 50uL of OS(=O)(=O)O
2080
2081 Chlorotrifluoromethane = 50uL of C(F)(F)(F)Cl
2082
2083 Dibutyl-phosphite = 50uL of CCCCOP(=O)(O)CCCC
2084
2085 Ammonium-dichromate = 50uL of [NH4+].[NH4+].[O-][Cr](=O)(=O)O[Cr](=O)(=O)[O-]
2086
2087
2088 a = Mix Sodium-Hydroxide with Chlorotrifluoromethane
2089 b = Mix a with Isopropyl-Alcohol
2090 c = Mix Allyl-Ethyl-Ether with Cycloheptatriene
2091 Mix Cycloheptatriene2 with Dibutyl-phosphite
2092

2093 7.8 Broad Spectrum Opiate

2094 Move 10uL of UrineSample to Anti-Morphine
2095 Move 10uL of UrineSample to Anti-Oxy
2096 Move 10uL of UrineSample to Anti-Fentanyl
2097 Move 10uL of UrineSample to Anti-Ciprofloxacin
2098 Move 10uL of UrineSample to Anti-Heroin
2099
2100 MorphineReading = Measure the fluorescence of Anti-Morphine for 5s
2101 OxyReading = Measure the fluorescence of Anti-Oxy for 5s
2102 FentanylReading = Measure the fluorescence of Anti-Fentanyl for 5s
2103 CiproReading = Measure the fluorescence of Anti-Ciprofloxacin for 5s
2104 HeroinReading = Measure the fluorescence of Anti-Heroin for 5s
2105
2106 Drain Anti-Morphine
2107

2108 Drain Anti-Oxy

2109 Drain Anti-Fentanyl

2110 Drain Anti-Ciprofloxacin

2111 Drain Anti-Heroin

2112

2113 7.9 Cipro ELISA

2114 Move 20uL of UrineSample to Ciprofloxacin-Enzyme

2115 Move 100uL of Ciprofloxacin-Conjugate to Ciprofloxacin-Enzyme

2116 Tap Ciprofloxacin-Enzyme for 60s

2117 Incubate Ciprofloxacin-Enzyme at 23 C for 60min

2118 Drain Ciprofloxacin-Enzyme

2119

2120 Repeat 5 times {

2121 Move 250uL of Distilled-Water to Ciprofloxacin-Enzyme

2122 Vortex Ciprofloxacin-Enzyme for 45s

2123 Drain Ciprofloxacin-Enzyme

2124 }

2125

2126 Move 50uL of TMB-Substrate to Ciprofloxacin-Enzyme

2127 Incubate Ciprofloxacin-Enzyme at 25 C for 30min

2128 Mix Ciprofloxacin-Enzyme with 100uL of Stop Reagent for 60s

2129

2130 Urine Reading = Measure the fluorescence of Ciprofloxacin-Enzyme for 5min

2131 Drain Ciprofloxacin-Enzyme

2132

2133 7.10 Diazepam ELISA

2134 Move 50uL of UrineSample to Diazepam-Enzyme

2135 Move 100uL of Diazepam-Antibody to Diazepam-Enzyme

2136 Tap Diazepam-Enzyme for 60s

2137 Incubate Diazepam-Enzyme at 23 C for 30min

2138 Drain Diazepam-Enzyme

2139

2140 Repeat 3 times {

2141 Move 250uL of Distilled-Water to Diazepam-Enzyme

2142 Vortex Diazepam-Enzyme for 45s

2143 Drain Diazepam-Enzyme

2144 }

2145

2146 Move 150uL of HRP-Conjugate to Diazepam-Enzyme

2147 Incubate Diazepam-Enzyme at 23 C for 15min

2148 Drain Diazepam-Enzyme

2149

2150 Repeat 3 times {

2151 Move 250uL of Distilled-Water to Diazepam-Enzyme

2152 Vortex Diazepam-Enzyme for 45s

2153 Drain Diazepam-Enzyme

2154 }

2155

2156

2157 Move 100uL of TMB-Substrate to Diazepam-Enzyme
2158 Incubate Diazepam-Enzyme at 23 C for 15min
2159 Mix Diazepam-Enzyme with 100uL of Stop Reagent for 60s
2160 Negative Reading = Measure the fluorescence of Diazepam-Enzyme for 30min
2161 Drain Diazepam-Enzyme
2162

2163 **7.11 547 Dilution**

2164 First-Dilute = Mix Substance A with Dilutant1
2165 Split First-Dilute into 2 parts:
2166 1 part to First-Dilute
2167 1 part to Waste-Droplet
2168 Drain Waste-Droplet
2169

2170 Second-Dilute = Mix First-Dilute with Dilutant2
2171 Split Second-Dilute into 2 parts:
2172 1 part to Second-Dilute
2173 1 part to Waste-Droplet
2174 Drain Waste-Droplet
2175

2176 Third Dilute = Mix Second-Dilute with Dilutant3
2177 Split Third-Dilute into 2 parts:
2178 1 part to Third-Dilute
2179 1 part to Waste-Droplet
2180 Drain Waste-Droplet
2181

2182 Fourth-Dilute = Mix Substance B with Substance C
2183 Split Fourth-Dilute into 2 parts:
2184 1 part to Fourth-Dilute
2185 1 part to Waste-Droplet
2186 Drain Waste-Droplet
2187

2188 Final-Dilute = Mix Third-Dilute with Fourth-Dilute
2189 Split Final-Dilute into 2 parts:
2190 1 part to Final-Dilute
2191 1 part to Waste-Droplet
2192 Drain Waste-Droplet
2193

2194 Save Final-Dilute
2195

2196 **7.12 Fentanyl ELISA**

2197
2198 # Step 1 Add Urine-Sample to test against
2199 Antigen1 = Move 20uL of Urine-Sample to Antigen
2200
2201 # Step 2 Add Enzyme Conjugate
2202 Antigen2 = Move 100uL of Fentanyl-Conjugate to Antigen1
2203
2204 # Step 3 Ensure fully mixed interactions
2205

```
2206 Tap Antigen2 for 60s
2207
2208 # Step 4
2209 Incubate Antigen2 at 23 C for 60min
2210
2211 # Step 5
2212 Drain Antigen2
2213 Repeat 6 times {
2214     Move 350uL of Distilled-Water to Antigen1
2215
2216     Vortex Antigen1 for 45s
2217
2218     Drain Antigen1
2219 }
2220 # Step 6
2221 Antigen1 = Move 100uL of TMB-Substrate to Antigen1
2222
2223 # Step 7
2224 Incubate Antigen1 at 23 C for 30min
2225
2226 # Step 8
2227 Mix Antigen1 with 100uL of Stop-Reagent for 60s
2228
2229 # Step 9
2230 Negative-Reading = Measure the fluorescence of Antigen1 for 30min
2231
2232 # End
2233 Drain Antigen1
```

7.13 Morphine ELISA

```
2235 Stationary Morphine-Enzyme
2236
2237 Move 20uL of Urine-Sample to Morphine-Enzyme
2238 Move 100uL of Morphine-Conjugate to Morphine-Enzyme
2239 Tap Morphine-Enzyme for 60s
2240 Incubate Morphine-Enzyme at 23 C for 60min
2241
2242 Repeat 6 times {
2243     Move 350uL of Distilled-Water to Morphine-Enzyme
2244     Vortex Morphine-Enzyme for 45s
2245 }
2246
2247 Move 100uL of TMB-Substrate to Morphine-Enzyme
2248 Incubate Morphine-Enzyme at 23 C for 30min
2249 Mix Morphine-Enzyme with 100uL of Stop Reagent for 60s
2250
2251 Urine Reading = Measure the fluorescence of Morphine-Enzyme for 30min
2252
2253
2254
```

7.14 Morphine ELISA with Control Samples

Step 2

Move 20uL of Negative-Standard to Antigen1

Move 20uL of Positive-Standard to Antigen2

Move 20uL of Diluted-Sample to Antigen3

Step 3

Move 100uL of Morphine-Conjugate to Antigen1

Move 100uL of Morphine-Conjugate to Antigen2

Move 100uL of Morphine-Conjugate to Antigen3

Step 4

Tap Antigen1 for 60s

Tap Antigen2 for 60s

Tap Antigen3 for 60s

Step 5

Incubate Antigen1 at 23 C for 60min

Incubate Antigen2 at 23 C for 60min

Incubate Antigen3 at 23 C for 60min

Step 6

Drain Antigen1

Drain Antigen2

Drain Antigen3

Repeat 6 times {

Move 350uL of Distilled-Water to Antigen1

Move 350uL of Distilled-Water to Antigen2

Move 350uL of Distilled-Water to Antigen3

Vortex Antigen1 for 45s

Vortex Antigen2 for 45s

Vortex Antigen3 for 45s

Drain Antigen1

Drain Antigen2

Drain Antigen3

}

Step 7 Dry...

Step 8

Move 100uL of TMB-Substrate to Antigen1

Move 100uL of TMB-Substrate to Antigen2

Move 100uL of TMB-Substrate to Antigen3

Step 9

2304 Incubate Antigen1 at 23 C for 30min
2305 Incubate Antigen2 at 23 C for 30min
2306 Incubate Antigen3 at 23 C for 30min
2307
2308 # Step 10
2309 Mix Antigen1 with 100uL of Stop Reagent for 60s
2310 Mix Antigen2 with 100uL of Stop Reagent for 60s
2311 Mix Antigen3 with 100uL of Stop Reagent for 60s
2312
2313 # Step 11
2314 Negative Reading = Measure the fluorescence of Antigen1 for 30min
2315 Positive Reading = Measure the fluorescence of Antigen2 for 30min
2316 Sample Reading = Measure the fluorescence of Antigen3 for 30min
2317
2318 # End
2319 Drain Antigen1
2320 Drain Antigen2
2321 Drain Antigen3
2322
2323 **7.15 Heroin ELISA**
2324 Move 20uL of Urine-Sample to Heroin-Enzyme
2325 Move 100uL of Heroin-Conjugate to Heroin-Enzyme
2326 Tap Heroin-Enzyme for 60s
2327 Incubate Heroin-Enzyme at 23 C for 60min
2328 Drain Heroin-Enzyme
2329
2330 Repeat 6 times {
2331 Move 350uL of Distilled-Water to Heroin-Enzyme
2332 Vortex Heroin-Enzyme for 45s
2333 Drain Heroin-Enzyme
2334 }
2335
2336 Move 100uL of TMB-Substrate to Heroin-Enzyme
2337 Incubate Heroin-Enzyme at 23 C for 30min
2338 Mix Heroin-Enzyme with 100uL of Stop Reagent for 60s
2339
2340 Urine-Reading = Measure the fluorescence of Heroin-Enzyme for 30min
2341 Drain Heroin-Enzyme
2342
2343 **7.16 Oxycodone ELISA**
2344 Move 20uL of Urine-Sample to Oxycodone-Enzyme
2345 Move 100uL of Oxycodone-Conjugate to Oxycodone-Enzyme
2346 Tap Oxycodone-Enzyme for 60s
2347 Incubate Oxycodone-Enzyme at 23 C for 60min
2348 Drain Oxycodone-Enzyme
2349
2350 Repeat 6 times {
2351 Move 350uL of Distilled-Water to Oxycodone-Enzyme
2352

2353 Vortex Oxycodone-Enzyme for 45s

2354 Drain Oxycodone-Enzyme

2355 }

2356

2357 Move 100uL of TMB-Substrate to Oxycodone-Enzyme

2358 Incubate Oxycodone-Enzyme at 23 C for 30min

2359 Mix Oxycodone-Enzyme with 100uL of Stop Reagent for 60s

2360

2361 Urine Reading = Measure the fluorescence of Oxycodone-Enzyme for 30min

2362 Drain Oxycodone-Enzyme

2363

2364 7.17 Image Probe Synthesis

2365 Mixture = Mix 10uL of Ion exchange beads with 10uL of Fluoride ions F⁻ for 30s

2366

2367 Heat Mixture at 100 C for 30s

2368 Heat Mixture at 120 C for 30s

2369 Heat Mixture at 135 C for 3min

2370

2371 Mixture = Mix Mixture with 10uL of MeCN solution for 30s

2372

2373 Incubate Mixture at 100 C for 30s

2374 Incubate Mixture at 120 C for 50s

2375

2376 Mixture = Mix Mixture with 10uL of HCl for 60s

2377 Heat Mixture at 60 C for 60s

2378

2379 7.18 Glucose Detection

2380 Result1 = Mix 10uL of Glucose with 10uL of Reagent for 10s

2381 Reading1 = Measure the fluorescence of Result1 for 30s

2382

2383 Result2 = Mix 10uL of Glucose with 20uL of Reagent for 10s

2384 Reading2 = Measure the fluorescence of Result2 for 30s

2385

2386 Result3 = Mix 10uL of Glucose with 40uL of Reagent for 10s

2387 Reading3 = Measure the fluorescence of Result3 for 30s

2388

2389 Result4 = Mix 10uL of Glucose with 80uL of Reagent for 10s

2390 Reading4 = Measure the fluorescence of Result4 for 30s

2391

2392 Result5 = Mix 10uL of Sample with 10uL of Reagent for 10s

2393 Reading5 = Measure the fluorescence of Result5 for 30s

2394

2395 7.19 PCR

2396 PCR-Mixture = 10uL of PCR

2397

2398 Heat PCR-Mixture at 95 C for 5s

2399

2400 Repeat 20 times {

2401

```

2402     Heat PCR-Mixture at 53 C for 15s
2403     Heat PCR-Mixture at 72 C for 10s
2404 }
2405 Perform Capillary Electrophoresis ( 5 cm at 236 V/cm) on PCR-Mixture Separate with Separati
2406
2407 Measure the fluorescence of PCR-Mixture for 3min
2408
2409 7.20 Neurotransmitter Sensing
2410 Mixture = Mix Sample with Reagent for 50s
2411 Perform Capillary Electrophoresis ( 9 cm at 223 V/cm) on Mixture Separate with electrophore
2412 Measure the fluorescence of Mixture for 10s
2413
2414 8 AQUACORE ASSAYS
2415 8.1 Glucose Detection
2416
2417 Input port ip1 ;Standard glucose
2418 Input port ip2 ;Reagent
2419 Input port ip3 ;Sample
2420 RESULT[1..5] ; dry array for final results
2421 glucose-detection {
2422   input s1, ip1
2423   input s2, ip2
2424   input s3, ip3
2425   move mixer1, s1, 1;5s
2426   move mixer1, s2, 1;5s
2427   mix mixer1, 10;10s
2428   move sensor1, mixer1;5s
2429   sense.OD sensor1, RESULT[1] ;30s
2430
2431   move mixer1, s1, 1 ;5s
2432   move mixer1, s2, 2;5s
2433   mix mixer1, 10;10s
2434   move sensor1, mixer1;5s
2435   sense.OD sensor1, RESULT[2] ;30s
2436
2437   move mixer1, s1, 1 ;5s
2438   move mixer1, s2, 4 ;5s
2439   mix mixer1, 10 ;10s
2440   move sensor1, mixer1 ;5s
2441   sense.OD sensor1, RESULT[3] ;30s
2442
2443   move mixer1, s1, 1 ;5s
2444   move mixer1, s2, 8 ;5s
2445   mix mixer1, 10 ;10s
2446   move sensor1, mixer1 ;5s
2447   sense.OD sensor1, RESULT[4] ;30s
2448
2449 <dry routine to get best line fit for RESULT[1..4]> move mixer1, s3, 1 ;5s
2450

```

```
2451 move mixer1, s2, 1 ;5s
2452 mix mixer1, 10 ;10s
2453 move sensor1, mixer1 ;5s
2454 sense.OD sensor1, RESULT[5] ;30s
2455
2456 <dry routine to get concentration from line given RESULT[5]>
2457 }
2458 Total time = 275s
2459 #reservoirs = 3
2460 ASLoC
```

2461 8.2 PCR

```
2462
2463 Input port ip1 ;PCR mixture
2464 Input port ip2 ;CE separation medium
2465
2466 RESULT[] ;dry array for final results
2467
2468 PCR {
2469   input s1, ip1
2470   input s2, ip2
2471
2472   move heater1, s1 ;5s
2473   dry-mov r1, 20
2474   dry-label loop:
2475     incubate heater1, 95, 5 ;6s
2476     incubate heater1, 53, 15 ;17s
2477     incubate heater1, 72, 10 ;12s
2478     dry-dec r1
2479     dry-bgt loop
2480
2481     move separator1.buf, s2 ;5s
2482     move separator1, heater1 ;5s
2483
2484     separate.CE separator1, 236, 5, 180
2485     sense.FL sensor1, RESULT ;180s
2486   }
2487   Total time = 895s
2488   #reservoirs = 2
2489   ASLoC
```

2491 8.3 Imaging Probe Synthesis

```
2492 Input port ip1 ;Ion exchange beads (in buffer)
2493 Input port ip2 ;Fluoride ions F-
2494 Input port ip3 ;MeCN solution
2495 Input port ip4 ;HCl
2496
2497 Imaging-Probe-Synthesis {
2498   input s1, ip1
```

```
2499
```

```

2500 input s2, ip2
2501 input s3, ip3
2502 input s4, ip4
2503
2504 move mixer1, s1
2505 move mixer1, s2
2506 mix mixer1, 30
2507
2508 move heater1, mixer1
2509 concentrate.EV heater1, 100, 30
2510 concentrate.EV heater1, 120, 30
2511 concentrate.EV heater1, 135, 180
2512
2513 move mixer1, s3
2514 move mixer1, heater1
2515 mix mixer1, 30
2516
2517 move heater1, mixer1
2518 incubate heater1, 100, 30
2519 incubate heater1, 120, 50
2520 move mixer1, heater1
2521
2522 move mixer1, s4
2523 mix mixer1, 60
2524 move heater1, mixer1
2525 concentrate.EV heater1, 60, 60
2526 }
2527 Total time = 548s #reservoirs = 4
2528 ASLoC area = 13mm x 9mm = 117 mm2
2529
2530 8.4 Neurotransmitter Sensing
2531 Input port ip1 ;Sample
2532 Input port ip2 ;Reagent(OPA)
2533 Input port ip3 ;Electrophoresis buffer
2534
2535 RESULT[] ;dry array for final results
2536
2537 neurotransmitter-sensing {
2538 input s1, ip1
2539 input s2, ip2
2540 input s3, ip3
2541
2542 move mixer1, s1 ;5s
2543 move mixer1, s2 ;5s
2544 mix mixer1, 50 ;50s
2545
2546 move separator1.buf, s3 ;5s
2547 move separator1, mixer1 ;5s
2548

```

```

2549
2550 separate.CE separator1, 223, 9, 22
2551 sense.FL sensor1, RESULT ;22s
2552
2553 }
2554 Total time = 92s #reservoirs = 3
2555 ASLoC area = 2.5cm x 1.5cm = 3.75cm2 + CE column
2556 FIGURE 3: Separations-based sensing of neurotransmitters (immunoassay, biochemistry): (a)
2557 tio
2558

```

2559 9 BIOCODER ASSAYS

2560 9.1 PCR

```

2561 void PCR(){
2562   BioSystem bioCoder;
2563
2564   Fluid *PCRMix = bioCoder.new_fluid("PCRMasterMix", Volume(MICRO_LITER,10));
2565   Fluid *SeperationMedium = bioCoder.new_fluid("Seperation Medium",Volume(MICRO_LITER,10));
2566   Container* tube = bioCoder.new_container(STERILE_MICROFUGE_TUBE2ML);
2567
2568   bioCoder.first_step();
2569   bioCoder.measure_fluid(PCRMix,tube);
2570
2571   bioCoder.next_step();
2572   bioCoder.incubate(tube,95,Time(5,SECS));
2573
2574   bioCoder.next_step();
2575   bioCoder.LOOP(20);
2576
2577   bioCoder.next_step();
2578   bioCoder.incubate(tube,53,Time(15,SECS));
2579
2580   bioCoder.next_step();
2581   bioCoder.incubate(tube,72,Time(10,SECS));
2582   bioCoder.END_LOOP();
2583
2584   bioCoder.next_step();
2585   bioCoder.ce_detect(tube,5,236,SeperationMedium);
2586
2587   bioCoder.next_step();
2588   bioCoder.measure_fluorescence(tube,Time(3,MINS));
2589
2590   bioCoder.next_step();
2591   bioCoder.end_protocol();
2592 }
2593
2594 9.2 Probabilistic PCR
2595 void ProbablisticPCR()
2596
2597

```

```
2598 {
2599   BioSystem bioCoder;
2600
2601   Fluid *PCRMix = bioCoder.new_fluid("PCRMasterMix", Volume(MICRO_LITER,10));
2602
2603   Container* tube = bioCoder.new_container(STERILE_MICROFUGE_TUBE2ML);
2604
2605   bioCoder.first_step();
2606   bioCoder.measure_fluid(PCRMix,tube);
2607
2608   for(int i = 0 ; i < initial; ++i) {
2609     bioCoder.next_step();
2610     bioCoder.store_for(tube,94,Time(SECS,45));
2611
2612     bioCoder.next_step();
2613     bioCoder.store_for(tube,65,Time(SECS,45));
2614   }
2615
2616   for(int i = initial; i <= Threshold; ++i) {
2617     std::cout <<i<<std::endl;
2618     bioCoder.next_step();
2619     bioCoder.store_for(tube,94,Time(SECS,45));
2620
2621     bioCoder.next_step();
2622     bioCoder.store_for(tube,65,Time(SECS,45));
2623
2624     bioCoder.next_step();
2625     bioCoder.measure_fluorescence(tube,Time(SECS,5),"DNASensor");
2626
2627     bioCoder.IF("DNASensor",GREATER_THAN, .85);
2628     for(int j = i; j < Total+(Threshold-i); ++j) {
2629       bioCoder.next_step();
2630       bioCoder.store_for(tube,94,Time(SECS,45));
2631
2632       bioCoder.next_step();
2633       bioCoder.store_for(tube,65,Time(SECS,45));
2634     }
2635
2636     bioCoder.next_step();
2637     bioCoder.drain(tube,"Amplified PCR");
2638     bioCoder.END_IF();
2639
2640   }
2641
2642   bioCoder.drain(tube,"waste");
2643   bioCoder.end_protocol();
2644
2645   bioCoder.PrintLeveledProtocol();
2646
```

```

2647 bioCoder.PrintTree();
2648 bioCoder.PrintTreeVisualization("ProbablisticPCR");
2649 }
2650
2651 9.3 PCR Droplet Replenish
2652 void PCRDropletReplacement()
2653 {
2654     int TotalThermo = 9;
2655     BioSystem bioCoder;
2656
2657     Fluid *PCRMix = bioCoder.new_fluid("PCRMasterMix", Volume(MICRO_LITER,10));
2658     Fluid *Template = bioCoder.new_fluid("Template", Volume(MICRO_LITER,10));
2659
2660     Container* tube = bioCoder.new_container(STERILE_MICROFUGE_TUBE2ML);
2661     //Container* tube2 = bioCoder.new_container(STERILE_MICROFUGE_TUBE2ML);
2662
2663     bioCoder.first_step();
2664     bioCoder.measure_fluid(PCRMix,tube);
2665
2666     bioCoder.next_step();
2667     bioCoder.vortex(tube,Time(SECS,1));
2668     bioCoder.measure_fluid(Template,tube);
2669
2670     bioCoder.next_step();
2671     bioCoder.vortex(tube, Time(SECS,1));
2672
2673     bioCoder.next_step();
2674     bioCoder.store_for(tube,95,Time(SECS,45));
2675
2676     bioCoder.next_step();
2677     bioCoder.LOOP(TotalThermo);
2678
2679     std::cout<<"Debug statement2"<<std::endl;
2680     bioCoder.next_step();
2681     bioCoder.store_for(tube,95,Time(SECS,20));
2682
2683     bioCoder.next_step();
2684     bioCoder.weigh(tube,"weightSensor");
2685
2686     bioCoder.next_step();
2687     bioCoder.IF("WieghtSensor",LESS_THAN, 3.57);
2688     bioCoder.next_step();
2689     bioCoder.measure_fluid(PCRMix, tube);
2690
2691     bioCoder.next_step();
2692     bioCoder.store_for(tube, 95,Time(SECS,45));
2693
2694     bioCoder.next_step();
2695

```

```

2696 bioCoder.vortex(tube, Time(SECS,1));
2697 bioCoder.END_IF();
2698
2699 bioCoder.next_step();
2700 bioCoder.store_for(tube,50,Time(SECS,30));
2701
2702 bioCoder.next_step();
2703 bioCoder.store_for(tube,68,Time(SECS,45));
2704 std::cout<<"Debug statement3"<<std::endl;
2705 bioCoder.END_LOOP();
2706 std::cout<<"Debug statement4"<<std::endl;
2707
2708
2709 bioCoder.next_step();
2710 bioCoder.store_for(tube,68,Time(MINS,5));
2711 std::cout<<"Debug statement5"<<std::endl;
2712 bioCoder.next_step();
2713 bioCoder.drain(tube,"PCR");
2714 std::cout<<"Debug statement6"<<std::endl;
2715 bioCoder.end_protocol();
2716
2717
2718 std::cout<<"Debug statemen7"<<std::endl;
2719 bioCoder.PrintLeveledProtocol();
2720 bioCoder.PrintTree();
2721 bioCoder.PrintTreeVisualization("PCRReplacement");
2722
2723 }

```

9.4 Glucose Detection

```

2726 void GlucoseDetection(){
2727 BioSystem bioCoder;
2728
2729 Fluid *Glucose = bioCoder.new_fluid("Ion exchange beads", Volume(MICRO_LITER,160));
2730 Fluid *Reagent = bioCoder.new_fluid("Fluoride ions",Volume(MICRO_LITER,50));
2731 Fluid *Sample = bioCoder.new_fluid("HCL",Volume(MICRO_LITER,10));
2732
2733 Container* tube = bioCoder.new_container(STERILE_MICROFUGE_TUBE2ML);
2734 Container* tube2 = bioCoder.new_container(STERILE_MICROFUGE_TUBE2ML);
2735 Container* tube3 = bioCoder.new_container(STERILE_MICROFUGE_TUBE2ML);
2736 Container* tube4 = bioCoder.new_container(STERILE_MICROFUGE_TUBE2ML);
2737 Container* tube5 = bioCoder.new_container(STERILE_MICROFUGE_TUBE2ML);
2738
2739 bioCoder.first_step();
2740 bioCoder.measure_fluid(Glucose,Volume(MICRO_LITER,10),tube);
2741 bioCoder.measure_fluid(Glucose,Volume(MICRO_LITER,10),tube2);
2742 bioCoder.measure_fluid(Glucose,Volume(MICRO_LITER,10),tube3);
2743 bioCoder.measure_fluid(Glucose,Volume(MICRO_LITER,10),tube4);
2744

```



```
2745 bioCoder.measure_fluid(Glucose,Volume(MICRO_LITER,10),tube5);
2746
2747 bioCoder.measure_fluid(Reagent,Volume(MICRO_LITER,10),tube);
2748 bioCoder.measure_fluid(Reagent,Volume(MICRO_LITER,20),tube2);
2749 bioCoder.measure_fluid(Reagent,Volume(MICRO_LITER,40),tube3);
2750 bioCoder.measure_fluid(Reagent,Volume(MICRO_LITER,80),tube4);
2751
2752 bioCoder.measure_fluid(Sample,Volume(MICRO_LITER,80),tube5);
2753
2754 bioCoder.next_step();
2755 bioCoder.measure_fluorescence(tube,Time(5,SECS));
2756 bioCoder.measure_fluorescence(tube2,Time(5,SECS));
2757 bioCoder.measure_fluorescence(tube3,Time(5,SECS));
2758 bioCoder.measure_fluorescence(tube4,Time(5,SECS));
2759 bioCoder.measure_fluorescence(tube5,Time(5,SECS));
2760
2761 bioCoder.next_step();
2762 bioCoder.end_protocol();
2763 }
```

2765 9.5 Image Probe Synthesis

```
2766 void ImageProbSynthesis(){
2767     BioSystem bioCoder;
2768
2769     Fluid *IonBeads = bioCoder.new_fluid("Ion exchange beads", Volume(MICRO_LITER,10));
2770     Fluid *Fluoride = bioCoder.new_fluid("Fluoride ions",Volume(MICRO_LITER,10));
2771     Fluid *HCL = bioCoder.new_fluid("HCL",Volume(MICRO_LITER,10));
2772     Fluid *MeCNSolution = bioCoder.new_fluid("MeCN solution",Volume(MICRO_LITER,10));
2773
2774     Container* tube = bioCoder.new_container(STERILE_MICROFUGE_TUBE2ML);
2775
2776     bioCoder.first_step();
2777     bioCoder.measure_fluid(IonBeads,tube);
2778     bioCoder.measure_fluid(Fluoride,tube);
2779
2780     bioCoder.next_step();
2781     bioCoder.vortex(tube,Time(30,SECS));
2782
2783     bioCoder.next_step();
2784     bioCoder.incubate(tube,100,Time(30,SECS));
2785
2786     bioCoder.next_step();
2787     bioCoder.incubate(tube,120,Time(30,SECS));
2788
2789     bioCoder.next_step();
2790     bioCoder.incubate(tube,135, Time(3,MINS));
2791
2792     bioCoder.next_step();
2793 }
```

```
2794 bioCoder.measure_fluid(MeCNSolution,tube);
2795
2796 bioCoder.next_step();
2797 bioCoder.vortex(tube,Time(30,SECS));
2798
2799 bioCoder.next_step();
2800 bioCoder.incubate(tube,100,Time(30,SECS));
2801
2802 bioCoder.next_step();
2803 bioCoder.incubate(tube,120,Time(50,SECS));
2804
2805 bioCoder.next_step();
2806 bioCoder.measure_fluid(HCL,tube);
2807
2808 bioCoder.next_step();
2809 bioCoder.vortex(tube,Time(60,SECS));
2810
2811 bioCoder.next_step();
2812 bioCoder.incubate(tube,60,Time(60,SECS));
2813
2814 bioCoder.next_step();
2815 bioCoder.end_protocol();
2816 }
```

2817 **9.6 Neurotransmitter Sensing**

```
2818 void neurotransmitterSensing(){
2819 BioSystem bioCoder;
2820
2821
2822 Fluid *Sample = bioCoder.new_fluid("Sample", Volume(MICRO_LITER,10));
2823 Fluid *Reagent = bioCoder.new_fluid("Reagent",Volume(MICRO_LITER,10));
2824 Fluid *SeperationMedium = bioCoder.new_fluid("Seperation Medium",Volume(MICRO_LITER,10));
2825
2826 Container* tube = bioCoder.new_container(STERILE_MICROFUGE_TUBE2ML);
2827
2828 bioCoder.first_step();
2829 bioCoder.measure_fluid(Sample,tube);
2830 bioCoder.measure_fluid(Reagent,tube);
2831
2832 bioCoder.next_step();
2833 bioCoder.vortex(tube,Time(50,SECS));
2834
2835 bioCoder.next_step();
2836 bioCoder.ce_detect(tube,9,223,SeperationMedium);
2837
2838 bioCoder.next_step();
2839 bioCoder.measure_fluorescence(tube,Time(10,SECS));
2840
2841 bioCoder.next_step();
2842
```

```

2843 bioCoder.end_protocol();
2844 }
2845
2846 10 ANTHA ASSAYS
2847 10.1 Glucose Detection
2848
2849 protocol Glucose_Detection
2850
2851
2852 import (
2853 "github.com/antha-lang/antha/antha/anthalib/wtype"
2854     // LHComponent type
2855 "github.com/antha-lang/antha/antha/anthalib/wutil"
2856 "github.com/antha-lang/antha/antha/anthalib/mixer"
2857     //sample function is imported from mixed
2858
2859 // Input parameters
2860 Parameters (
2861 Reagent_volume1 Volume // 10ul
2862 Reagent_volume2 Volume // 20ul
2863 Reagent_volume3 Volume // 40ul
2864 Reagent_volume4 Volume // 80ul
2865 Glucose_volume Volume // 10ul
2866 Sample_volume Volume // 10ul
2867 )
2868
2869
2870 // Data which is returned from this protocol, and data types
2871 Data (
2872
2873 )
2874
2875 // Physical Inputs to this protocol with types
2876 Inputs (
2877 Glucose *wtype.LHComponent
2878 Reagent *wtype.LHComponent
2879 Sample *wtype.LHComponent
2880 )
2881
2882
2883
2884 // Physical outputs from this protocol with types
2885 Outputs (
2886 Result1 *wtype.LHComponent
2887 Result2 *wtype.LHComponent
2888 Result3 *wtype.LHComponent
2889 Result4 *wtype.LHComponent
2890 Result5 *wtype.LHComponent
2891

```

```

2892 )
2893
2894
2895 Requirements {
2896 }
2897
2898 // Conditions to run on startup
2899 Setup {
2900
2901 }
2902
2903 // The core process for this protocol, with the steps to be performed
2904 // for every input
2905 Steps {
2906
2907   glucose := mixer.Sample(Glucose, Glucose_volume)
2908   reagent := mixer.Sample(Reagent, Reagent_volume1)
2909
2910   Result1 = mixer.Mix(glucose, reagent) // cannot specify duration of mixture
2911   //cannot measure fluorescence to get a reading
2912
2913   glucose := mixer.Sample(Glucose, Glucose_volume) //get new sample of Glucose
2914   reagent := mixer.Sample(Reagent, Reagent_volume2) //get new sample of Reagent
2915
2916   Result2 = mixer.Mix(glucose, reagent) // mix new samples
2917   // would measure fluorescence here
2918
2919   glucose := mixer.Sample(Glucose, Glucose_volume) //get new sample of Glucose
2920   reagent := mixer.Sample(Reagent, Reagent_volume3) //get new sample of Reagent
2921
2922   Result3 = mixer.Mix(glucose, reagent) // mix new samples
2923   // would measure fluorescence here
2924
2925   glucose := mixer.Sample(Glucose, Glucose_volume) //get new sample of Glucose
2926   reagent := mixer.Sample(Reagent, Reagent_volume4) //get new sample of Reagent
2927
2928   Result4 = mixer.Mix(glucose, reagent) // mix new samples
2929   // would measure fluorescence here
2930
2931   sample := mixer.Sample(Sample, Sample_volume) //get new sample of Sample
2932   reagent := mixer.Sample(Reagent, Reagent_volume1) //get new sample of Reagent
2933
2934   Result5 = mixer.Mix(sample, reagent) // mix new samples
2935   // would measure fluorescence here
2936
2937 }
2938
2939 Analysis {
2940

```

```
2941
2942 }
2943
2944 Validation {
2945 }
2946
2947 10.2 Imaging Probe Synthesis
2948 protocol Imaging_Probe_Synthesis
2949
2950
2951 import (
2952 "github.com/antha-lang/antha/antha/anthalib/wtype" // LHComponent type
2953 "github.com/antha-lang/antha/antha/anthalib/wutil"
2954 "github.com/antha-lang/antha/antha/anthalib/mixer"
2955     //sample function is imported from mixed
2956
2957 // Input parameters
2958 Parameters (
2959 Ion_volume Volume
2960 Flouride_volume Volume
2961 MeCN_volume Volume
2962 HCl_volume Volume
2963 )
2964
2965
2966 // Data which is returned from this protocol, and data types
2967 Data (
2968
2969 )
2970
2971 // Physical Inputs to this protocol with types
2972 Inputs (
2973 Ion_exchange_beads *wtype.LHComponent
2974 Fluoride_ions_F- *wtype.LHComponent
2975 MeCN *wtype.LHComponent
2976 HCl *wtype.LHComponent
2977 )
2978
2979
2980 // Physical outputs from this protocol with types
2981 Outputs (
2982 Mixture *wtype.LHComponent
2983 )
2984
2985 Requirements {
2986 }
2987
2988 // Conditions to run on startup
2989
```

```

2990 Setup {
2991
2992 }
2993
2994 // The core process for this protocol, with the steps to be performed
2995 // for every input
2996 Steps {
2997
2998   ibv := mixer.Sample(Ion_exchange_beads, Ion_volume)
2999   fif := mixer.Sample(Flourid_ions_F-, Flouride_volume)
3000
3001   mixture := mixer.Mix(ibv, fif) //cannot specify time for mixture
3002
3003   mixture = Incubate(mixture, 100, 30, false) //heat 100 for 30s
3004   mixture = Incubate(mixture, 120, 30, false) //heat 120 for 30s
3005   mixture = Incubate(mixture, 135, 180, false) //heat 135 for 3 minutes
3006
3007   mecn := mixer.Sample(MeCN, MeCN_volume)
3008
3009   mixture = mixer.Mix(mixture, mecn)
3010
3011   mixture = Incubate(mixture, 100, 30, false) //incubate 100 for 30s
3012   mixture = Incubate(mixture, 120, 50, false) //incubate 120 for 50s
3013
3014   hcl := mixer.Sample(HCl, HCl_volume)
3015
3016   mixture = mixer.Mix(mixture, hcl)
3017
3018   mixture = Incubate(mixture, 60, 60, false) // heat 60 for 60s
3019
3020   Mixture = mixture
3021 }
3022
3023 Analysis {
3024
3025 }
3026
3027 Validation {
3028 }
3029
3030 10.3 PCR
3031 protocol PCR
3032
3033 import (
3034   "github.com/antha-lang/antha/antha/anthalib/wtype"
3035   "github.com/antha-lang/antha/antha/anthalib/mixer"
3036   "github.com/antha-lang/antha/antha/AnthaStandardLibrary/Packages/enzymes"
3037   "fmt"
3038

```

```
3039 "github.com/antha-lang/antha/antha/AnthaStandardLibrary/Packages/text"
3040 "github.com/antha-lang/antha/antha/AnthaStandardLibrary/Packages/search"
3041 "github.com/antha-lang/antha/antha/AnthaStandardLibrary/Packages/sequences"
3042 )
3043
3044 /*type Polymerase struct {
3045     wtype.LHComponent
3046     Rate_BPpers float64
3047     Fidelity_errorate float64 // could dictate how many colonies are checked in validation!
3048     Extensiontemp Temperature
3049     Hotstart bool
3050     StockConcentration Concentration // this is normally in U?
3051     TargetConcentration Concentration
3052     // this is also a glycerol solution rather than a watersolution!
3053 }
3054 */
3055
3056
3057 // Input parameters for this protocol (data)
3058 Parameters (
3059     // PCRprep parameters:
3060     ReactionVolume Volume
3061     FwdPrimerConc Concentration
3062     RevPrimerConc Concentration
3063     Additiveconc Concentration
3064     TargetpolymeraseConcentration Concentration
3065     Templatevolume Volume
3066     Dntpconc Concentration
3067     /*
3068     // let's be ambitious and try this as part of type polymerase Polymeraseconc Volume
3069
3070     //Templatetype string // e.g. colony, genomic, pure plasmid...
3071     //will effect efficiency. We could get more sophisticated here later on...
3072     //FullTemplatesequence string
3073     // better to use Sid's type system here after proof of concept
3074     //FullTemplatelength int
3075     // clearly could be calculated from the sequence...
3076     //Sid will have a method to do this already so check!
3077     //TargetTemplatesequence string
3078     // better to use Sid's type system here after proof of concept
3079     //TargetTemplatelengthinBP int
3080     */
3081     // Reaction parameters: (could be a entered as thermocycle parameters type possibly?)
3082     Numberofcycles int
3083     InitDenaturationtime Time
3084     Denaturationtime Time
3085     //Denaturationtemp Temperature
3086     Annealingtime Time
3087
```

```

3088 AnnealingTemp Temperature // Should be calculated from primer and template binding
3089 Extensiontime Time // should be calculated from template length and polymerase rate
3090 Extensiontemp Temperature
3091 Finalexextensiontime Time
3092 Targetsequence string
3093 FwdPrimerSeq string
3094 RevPrimerSeq string
3095 )
3096
3097 // Data which is returned from this protocol, and data types
3098 Data (
3099   FwdPrimerSites []search.Thingfound
3100   RevPrimerSites []search.Thingfound
3101 )
3102
3103
3104 // Physical Inputs to this protocol with types
3105 Inputs (
3106   FwdPrimer *wtype.LHComponent
3107   RevPrimer *wtype.LHComponent
3108   DNTPS *wtype.LHComponent
3109   PCRPolymerase *wtype.LHComponent
3110   Buffer *wtype.LHComponent
3111   Template *wtype.LHComponent
3112   Additives []*wtype.LHComponent // e.g. DMSO
3113   OutPlate *wtype.LHPlate
3114 )
3115
3116 // Physical outputs from this protocol with types
3117 Outputs (
3118   Reaction *wtype.LHComponent
3119 )
3120
3121 Requirements {
3122 }
3123
3124 // Conditions to run on startup
3125 Setup {
3126 }
3127
3128 // The core process for this protocol, with the steps to be performed
3129 // for every input
3130 Steps {
3131
3132   FwdPrimerSites = sequences.FindSeqsinSeqs(Targetsequence, []string{FwdPrimerSeq})
3133
3134   RevPrimerSites = sequences.FindSeqsinSeqs(Targetsequence, []string{RevPrimerSeq})
3135
3136

```



```

3137
3138 if len(FwdPrimerSites)==0 || len(RevPrimerSites)==0{
3139
3140 errordescription := fmt.Sprintf(
3141 text.Print("FwdPrimerSitesfound:", fmt.Sprintf(FwdPrimerSites)),
3142 text.Print("RevPrimerSitesfound:", fmt.Sprintf(RevPrimerSites)),
3143 )
3144
3145 Errorf(errordescription)
3146 }
3147
3148
3149
3150 // Mix components
3151 samples := make([]*wtype.LHComponent, 0)
3152 bufferSample := mixer.SampleForTotalVolume(Buffer, ReactionVolume)
3153 samples = append(samples, bufferSample)
3154 templateSample := mixer.Sample(Template, Templatevolume)
3155 samples = append(samples, templateSample)
3156 dntpSample := mixer.SampleForConcentration(DNTPS, DNTPconc)
3157 samples = append(samples, dntpSample)
3158 FwdPrimerSample := mixer.SampleForConcentration(FwdPrimer, FwdPrimerConc)
3159 samples = append(samples, FwdPrimerSample)
3160 RevPrimerSample := mixer.SampleForConcentration(RevPrimer, RevPrimerConc)
3161 samples = append(samples, RevPrimerSample)
3162
3163 for _, additive := range Additives {
3164 additiveSample := mixer.SampleForConcentration(additive, Additiveconc)
3165 samples = append(samples, additiveSample)
3166 }
3167
3168 polySample := mixer.SampleForConcentration(PCRPolymerase, TargetpolymeraseConcentration)
3169 samples = append(samples, polySample)
3170 reaction := MixInto(OutPlate, "", samples...)
3171
3172 // thermocycle parameters called from enzyme lookup:
3173
3174 polymerase := PCRPolymerase.CName
3175
3176 extensionTemp := enzymes.DNAPolymerasetemps[polymerase]["extensiontemp"]
3177 meltingTemp := enzymes.DNAPolymerasetemps[polymerase]["meltingtemp"]
3178
3179
3180 // initial Denaturation
3181
3182 r1 := Incubate(reaction, meltingTemp, InitDenaturationtime, false)
3183
3184 for i:=0; i < Numberofcycles; i++ {
3185

```

```

3186
3187 // Denature
3188
3189 r1 = Incubate(r1, meltingTemp, Denaturationtime, false)
3190
3191 // Anneal
3192 r1 = Incubate(r1, AnnealingTemp, Annealingtime, false)
3193
3194 //extensiontime := TargetTemplatelengthinBP/PCRPolymerase.RateBPpers
3195 // we'll get type issues here so leave it out for now
3196
3197 // Extend
3198 r1 = Incubate(r1, extensionTemp, Extensiontime, false)
3199
3200 }
3201 // Final Extension
3202 r1 = Incubate(r1, extensionTemp, Finalexextensiontime, false)
3203
3204
3205 // all done
3206 Reaction = r1
3207 }
3208
3209 // Run after controls and a steps block are completed to
3210 // post process any data and provide downstream results
3211 Analysis {
3212 }
3213
3214 // A block of tests to perform to validate that the sample was processed correctly
3215 // Optionally, destructive tests can be performed to validate results on a
3216 // dipstick basis
3217 Validation {
3218 }
3219
3220
3221 10.4 Neurotransmitter Sensing
3222 protocol Neurotransmitter_Sensing
3223
3224
3225 import (
3226 "github.com/antha-lang/antha/antha/anthalib/wtype" // LHComponent type
3227 "github.com/antha-lang/antha/antha/anthalib/wutil"
3228 "github.com/antha-lang/antha/antha/anthalib/mixer" //sample function is imported from mixer
3229
3230 // Input parameters
3231 Parameters (
3232 Sample_volume //no val specified
3233 Reagent_volume //no val specified
3234

```

```
3235 electrophoresis_buffer_volume //no val specified
3236 )
3237
3238
3239 // Data which is returned from this protocol, and data types
3240 Data (
3241
3242 )
3243
3244 // Physical Inputs to this protocol with types
3245 Inputs (
3246 Sample *wtype.LHComponent
3247 Reagent *wtype.LHComponent
3248 electrophoresis_buffer *wtype.LHComponent
3249 )
3250
3251
3252
3253 // Physical outputs from this protocol with types
3254 Outputs (
3255 Mixture *wtype.LHComponent
3256 )
3257
3258
3259 Requirements {
3260 }
3261
3262 // Conditions to run on startup
3263 Setup {
3264
3265 }
3266
3267 // The core process for this protocol, with the steps to be performed
3268 // for every input
3269 Steps {
3270
3271 sample := mixer.Sample(Sample, Sample_volume)
3272 reagent := mixer.Sample(Reagent, Reagent_volume)
3273
3274 Mixture = mixer.Mix(sample, reagent) // cannot specify duration of mixture
3275 //cannot perform capillary electrophoresis
3276 //cannot measure fluorescence to get a reading
3277 }
3278
3279 Analysis {
3280
3281 }
3282
3283
```

Table 2. Experimental tests that validate ChemType; parentheses denote reactive group numbers assigned to chemical names. The AIHA and Safety Zone tests are real-world documented incidents that ChemType could have prevented. I & C respectively denote Incompatible errors (known to be dangerous) and Caution errors (likely to be dangerous) based on the EPA/NOAA reactive groups. This table includes both synthetic and real-world results.

Experiment Name	Outcome	Experiment Description
AIHA 1 [2]	FAIL - I	Mix Nitric Acid (2) and Tetrachloroethylene (17,28)
AIHA 2 [2]	FAIL - I	Mix Nitric Acid (2) and Methanol (4)
AIHA 3 [2]	FAIL - I	Mix Potassium Hydride (35, 21) and Diaminopropane (7)
Mustard Gas [1]	FAIL - I	Mix Calcium Hypo (1) and Dichlor (17)
Safety Zone [12]	FAIL - I	Mix Hydrogen Peroxide (44) and Sulfuric Acid (2), then mix Acetone (19)
Sytn - Fail 1	FAIL - I	Mix Nitric Acid (2) and Hydrochloric Acid (1)
Sytn - Fail 2	FAIL - C	Mix Hydrochloric Acid (1) and Isopropanol (4)
Sytn - Fail 3	FAIL - C	Mix Benzoin (4, 16, 19) and Urea (6), then mix formaldehyde (5, 76)
Sytn - Pass 1	PASS	Mix Benzene (16) and Urea (6)
Sytn - Pass 2	PASS	Mix Benzene (16) and Urea (6), then mix Tetraethoxysilane (58)

Validation {
}

11 CHEMTYPE RESULTS

11.1 Efficacy of ChemType

11.2 ChemType Runtime Evaluation

Table 3. Execution time and number of constraints for synthetic and real-world bioassays.

Benchmark	Compilation Time (sec)	Constraint Solving Time (sec)	Gathered constraints	Con-
AIHA 1 [2]	0.012	0.936	70	
AIHA 2 [2]	0.012	1.648	68	
AIHA 3 [2]	0.014	1.214	17	
Broad Spectrum Opi- ate [4, 24, 35]	0.011	0.887	11	
Ciprofloxacin [24]	0.023	1.722	14	
Diazepam [20]	0.024	1.007	14	
Dilution [20]	0.014	0.892	9	
Fentanyl [35]	0.018	0.900	13	
Full Morphine [20]	0.048	4.188	19	
Glucose Detection [3]	0.012	1.633	14	
Heroin [20]	0.020	1.553	13	
Image Probe Synthe- sis [3]	0.015	2.181	13	
Morphine [20]	0.018	1.026	13	
Mustard Gas [1]	0.015	1.433	83	
Neurotransmitter Sensing [3]	0.029	3.845	6	
Oxycodone [4]	0.026	0.959	13	
PCR [3]	0.032	3.534	8	
Safety Zone [12]	0.013	1.341	76	
Synthetic 1	0.031	1.507	196	
Synthetic 2	0.014	1.505	274	
Synthetic 3	0.012	0.944	95	
Synthetic 4	0.014	0.623	1	
Synthetic 5	0.026	1.350	10	
Synthetic 6	0.013	0.855	8	

REFERENCES

[1] 2016. Swimming pool chemical incident. <http://bit.ly/2gghGZL>. (2016). accessed: 2016-11-01.

[2] American Industrial Hygiene Association. 2016. <http://bit.ly/2eZtf1m>. (2016). Accessed: 2016-11-08.

[3] Ahmed M. Amin, Mithuna Thottethodi, T. N. Vijaykumar, Steven Wereley, and Stephen C. Jacobson. 2007. Aquacore: a programmable architecture for microfluidics. In *34th International Symposium on Computer Architecture (ISCA 2007)*, June 9-13, 2007, San Diego, California, USA. 254–265. <https://doi.org/10.1145/1250662.1250694>

[4] Ronald C Backer, Joseph R Monforte, and Alphonse Poklis. 2005. Evaluation of the DRI® oxycodone immunoassay for the detection of oxycodone in urine. *Journal of analytical toxicology* 29, 7 (2005), 675–677.

[5] E. Barsoukov and J.R. Macdonald. 2005. Impedence Spectroscopy: Theory, Experiments, and Applications. *Wiley-Interscience* (2005).

[6] Amar S. Basu. 2013. Droplet morphometry and velocimetry (DMV): a video processing software for time-resolved, label-free tracking of droplet parameters. *Lab Chip* 13 (2013), 1892–1901. Issue 10. <https://doi.org/10.1039/C3LC50074H>

[7] Karl-Friedrich Böhringer. 2006. Modeling and Controlling Parallel Tasks in Droplet-Based Microfluidic Systems. *IEEE Trans. on CAD of Integrated Circuits and Systems* 25, 2 (2006), 334–344. <https://doi.org/10.1109/TCAD.2005.855958>

[8] Ying-Han Chen, Chung-Lun Hsu, Li-Chen Tsai, Tsung-Wei Huang, and Tsung-Yi Ho. 2013. A Reliability-Oriented Placement Algorithm for Reconfigurable Digital Microfluidic Biochips Using 3-D Deferred Decision Making Technique. *IEEE Trans. on CAD of Integrated Circuits and Systems* 32, 8 (2013), 1151–1162. <https://doi.org/10.1109/TCAD.2013.2249558>

- [9] Minsik Cho and David Z. Pan. 2008. A High-Performance Droplet Routing Algorithm for Digital Microfluidic Biochips. *IEEE Trans. on CAD of Integrated Circuits and Systems* 27, 10 (2008), 1714–1724. <https://doi.org/10.1109/TCAD.2008.2003282>
- [10] P. Cooreman, R. Thoelen, J. Manca, M. vandeVen, V. Vermeeren, L. Michiels, M. Ameloot, and P. Wagner. 2005. Impedimetric immunosensors based on the conjugated polymer PPV. *Biosens. Bioelectron.* 20 (2005), 2151–2156. Issue 10.
- [11] Jie Ding, Krishnendu Chakrabarty, and Richard B. Fair. 2001. Scheduling of microfluidic operations for reconfigurable two-dimensional electrowetting arrays. *IEEE Trans. on CAD of Integrated Circuits and Systems* 20, 12 (2001), 1463–1468. <https://doi.org/10.1109/43.969439>
- [12] DA Dobbs, RG Bergman, and KH Theopold. 1990. Piranha solution explosion. (1990), 2–2 pages.
- [13] Ryan Fobel, Christian Fobel, and Aaron R. Wheeler. 2013. DropBot: An open-source digital microfluidic control system with precise control of electrostatic driving force and instantaneous drop velocity measurement. *Applied Physics Letters* 102, 19, Article 193513 (2013). <https://doi.org/10.1063/1.4807118>
- [14] Georges G. E. Gielen (Ed.). 2006. *Proceedings of the Conference on Design, Automation and Test in Europe, DATE 2006, Munich, Germany, March 6-10, 2006*. European Design and Automation Association, Leuven, Belgium. <http://ieeexplore.ieee.org/xpl/mostRecentIssue.jsp?punumber=11014>
- [15] J. Gong and C.J. Kim. 2008. Direct-referencing two-dimensional-array digital microfluidics using multilayer printed circuit board. *J. Microelectromech. Syst.* 17 (2008), 257–264. Issue 2.
- [16] Daniel Grissom and Philip Brisk. 2012. Path scheduling on digital microfluidic biochips. In *The 49th Annual Design Automation Conference 2012, DAC '12, San Francisco, CA, USA, June 3-7, 2012*, Patrick Groeneveld, Donatella Sciuto, and Soha Hassoun (Eds.). ACM, 26–35. <https://doi.org/10.1145/2228360.2228367>
- [17] Daniel Grissom and Philip Brisk. 2014. Fast Online Synthesis of Digital Microfluidic Biochips. *IEEE Trans. on CAD of Integrated Circuits and Systems* 33, 3 (2014), 356–369. <https://doi.org/10.1109/TCAD.2013.2290582>
- [18] Daniel Grissom, Christopher Curtis, and Philip Brisk. 2014. Interpreting Assays with Control Flow on Digital Microfluidic Biochips. *ACM Journal on Emerging Technologies (JETC) in Computing Systems* 10 (April 2014). Issue 3.
- [19] B. Hadwen, G. R. Broder, D. Morganti, A. Jacobs, C. Brown, J. R. Hector, Y. Kubota, and H. Morgan. 2012. Programmable large area digital microfluidic array with integrated droplet sensing for bioassays. *Lab Chip* 12, 18 (Sep 2012), 3305–3313.
- [20] Peter V Hornbeck. 1991. Enzyme-linked immunosorbent assays. *Current protocols in immunology* (1991), 2–1.
- [21] Kai Hu, Bang-Ning Hsu, Andrew Madison, Krishnendu Chakrabarty, and Richard B. Fair. 2013. Fault detection, real-time error recovery, and experimental demonstration for digital microfluidic biochips. In *Design, Automation and Test in Europe, DATE 13, Grenoble, France, March 18-22, 2013*, Enrico Macii (Ed.). EDA Consortium San Jose, CA, USA / ACM DL, 559–564. <https://doi.org/10.7873/DATE.2013.124>
- [22] Tsung-Wei Huang and Tsung-Yi Ho. 2009. A fast routability- and performance-driven droplet routing algorithm for digital microfluidic biochips. In *27th International Conference on Computer Design, ICCD 2009, Lake Tahoe, CA, USA, October 4-7, 2009*. 445–450. <https://doi.org/10.1109/ICCD.2009.5413119>
- [23] Tsung-Wei Huang, Chun-Hsien Lin, and Tsung-Yi Ho. 2010. A Contamination Aware Droplet Routing Algorithm for the Synthesis of Digital Microfluidic Biochips. *IEEE Trans. on CAD of Integrated Circuits and Systems* 29, 11 (2010), 1682–1695. <https://doi.org/10.1109/TCAD.2010.2062770>
- [24] Yousheng Jiang, Xuanyun Huang, Kun Hu, Wenjuan Yu, Xianle Yang, and Liquan Lv. 2011. Production and characterization of monoclonal antibodies against small haptens-ciprofloxacin. *African Journal of Biotechnology* 10, 65 (2011), 14342–14347.
- [25] Oliver Keszocze, Robert Wille, Krishnendu Chakrabarty, and Rolf Drechsler. 2015. A General and Exact Routing Methodology for Digital Microfluidic Biochips. In *Proceedings of the IEEE/ACM International Conference on Computer-Aided Design, ICCAD 2015, Austin, TX, USA, November 2-6, 2015*, Diana Marculescu and Frank Liu (Eds.). IEEE, 874–881. <https://doi.org/10.1109/ICCAD.2015.7372663>
- [26] Oliver Keszocze, Robert Wille, and Rolf Drechsler. 2014. Exact routing for digital microfluidic biochips with temporary blockages. In *The IEEE/ACM International Conference on Computer-Aided Design, ICCAD 2014, San Jose, CA, USA, November 3-6, 2014*, Yao-Wen Chang (Ed.). IEEE, 405–410. <https://doi.org/10.1109/ICCAD.2014.7001383>
- [27] Y. Li, H. Li, and R.J. Baker. 2014. Volume and concentration identification by using an electrowetting on dielectric device. *IEEE DCAS* (2014), 1–4.
- [28] C. Liao and S. Hu. 2011. Multiscale variation-aware techniques for high-performance digital microfluidic lab-on-a-chip component placement. *IEEE Trans Nanobioscience* 10, 1 (Mar 2011), 51–58.
- [29] G. Lippmann. 1875. *Relations entre les phénomènes électriques et capillaires*. Gauthier-Villars. <https://books.google.ch/books?id=IkqVNwAACAAJ>
- [30] Chia-Hung Liu, Kuang-Cheng Liu, and Juinn-Dar Huang. 2013. Latency-optimization synthesis with module selection for digital microfluidic biochips. In *2013 IEEE International SOC Conference, Erlangen, Germany, September 4-6, 2013*, Norbert Schuhmann, Kaijian Shi, and Nagi Naganathan (Eds.). IEEE, 159–164. <https://doi.org/10.1109/SOCC.2013>

- 6749681
- [31] L. Luan, R.D. Evans, N.M. Jokerst, and R.B. Fair. 2008. Integrated optical sensor in a digital microfluidic platform. *IEEE Sensors* 8 (2008), 628–635. Issue 5.
 - [32] Yan Luo, Bhargab B. Bhattacharya, Tsung-Yi Ho, and Krishnendu Chakrabarty. 2015. Design and Optimization of a Cyberphysical Digital-Microfluidic Biochip for the Polymerase Chain Reaction. *IEEE Trans. on CAD of Integrated Circuits and Systems* 34, 1 (2015), 29–42. <https://doi.org/10.1109/TCAD.2014.2363396>
 - [33] Elena Maftai, Paul Pop, and Jan Madsen. 2010. Tabu search-based synthesis of digital microfluidic biochips with dynamically reconfigurable non-rectangular devices. *Design Autom. for Emb. Sys.* 14, 3 (2010), 287–307. <https://doi.org/10.1007/s10617-010-9059-x>
 - [34] Elena Maftai, Paul Pop, and Jan Madsen. 2013. Module-Based Synthesis of Digital Microfluidic Biochips with Droplet-Aware Operation Execution. *JETC* 9, 1 (2013), 2. <https://doi.org/10.1145/2422094.2422096>
 - [35] Chi-Liang Mao, Keith D Zientek, Patrick T Colahan, Mei-Yueh Kuo, Chi-Ho Liu, Kuo-Ming Lee, and Chi-Chung Chou. 2006. Development of an enzyme-linked immunosorbent assay for fentanyl and applications of fentanyl antibody-coated nanoparticles for sample preparation. *Journal of pharmaceutical and biomedical analysis* 41, 4 (2006), 1332–1341.
 - [36] H. Moon, S.K. Cho, and C.J. Garrell, R.and Kim. 2002. Low voltage electrowetting-on-dielectric. *J. Appl. Phys.* 92 (2002), 4080–4087. Issue 7.
 - [37] Frieder Mugele and Jeanchristophe Baret. 2005. Electrowetting: from basics to applications. *Journal of Physics: Condensed Matter* 17 (2005), R705–R774.
 - [38] Miguel Angel Murran and Homayoun Najjaran. 2012. Capacitance-based droplet position estimator for digital microfluidic devices. *Lab Chip* 12 (2012), 2053–2059. Issue 11. <https://doi.org/10.1039/C2LC21241B>
 - [39] J. H. Noh, J. Noh, E. Kreit, J. Heikenfeld, and P. D. Rack. 2012. Toward active-matrix lab-on-a-chip: programmable electrofluidic control enabled by arrayed oxide thin film transistors. *Lab Chip* 12, 2 (Jan 2012), 353–360.
 - [40] Kenneth O’Neal, Daniel Grissom, and Philip Brisk. 2012. Force-Directed List Scheduling for Digital Microfluidic Biochips. In *20th IEEE/IFIP International Conference on VLSI and System-on-Chip, VLSI-SoC 2012, Santa Cruz, CA, USA, October 7-10, 2012*, Srinivas Katkoori, Matthew R. Guthaus, Ayse Kivildim Coskun, Andreas Burg, and Ricardo Reis (Eds.). IEEE, 7–11. <https://doi.org/10.1109/VLSI-SoC.2012.6378997>
 - [41] MG Pollack, AD Shenderov, and RB Fair. 2002. Electrowetting-based actuation of droplets for integrated microfluidics. *Lab on a Chip* 2, 2 (2002), 96–101.
 - [42] Andrew J. Ricketts, Kevin M. Irick, Narayanan Vijaykrishnan, and Mary Jane Irwin. 2006. Priority scheduling in digital microfluidics-based biochips, See [14], 329–334. <https://doi.org/10.1109/DATE.2006.244178>
 - [43] Pranab Roy, Hafizur Rahaman, and Parthasarathi Dasgupta. 2010. A novel droplet routing algorithm for digital microfluidic biochips. In *Proceedings of the 20th ACM Great Lakes Symposium on VLSI 2009, Providence, Rhode Island, USA, May 16-18 2010*, R. Iris Bahar, Fabrizio Lombardi, David Atienza, and Erik Brunvand (Eds.). ACM, 441–446. <https://doi.org/10.1145/1785481.1785583>
 - [44] Pranab Roy, Hafizur Rahaman, and Parthasarathi Dasgupta. 2012. Two-level clustering-based techniques for intelligent droplet routing in digital microfluidic biochips. *Integration* 45, 3 (2012), 316–330. <https://doi.org/10.1016/j.vlsi.2011.11.006>
 - [45] S. C. Shih, I. Barbulovic-Nad, X. Yang, R. Fobel, and A. R. Wheeler. 2013. Digital microfluidics with impedance sensing for integrated cell culture and analysis. *Biosens Bioelectron* 42 (Apr 2013), 314–320.
 - [46] Steve C. C. Shih, Ryan Fobel, Paresh Kumar, and Aaron R. Wheeler. 2011. A feedback control system for high-fidelity digital microfluidics. *Lab Chip* 11 (2011), 535–540. Issue 3. <https://doi.org/10.1039/C0LC00223B>
 - [47] Yong Jun Shin and Jeong Bong Lee. 2010. Machine vision for digital microfluidics. *Review of Scientific Instruments* 81, 1 (2 2010). <https://doi.org/10.1063/1.3274673>
 - [48] Fei Su and Krishnendu Chakrabarty. 2006. Module placement for fault-tolerant microfluidics-based biochips. *ACM Trans. Design Autom. Electr. Syst.* 11, 3 (2006), 682–710. <https://doi.org/10.1145/1142980.1142987>
 - [49] Fei Su and Krishnendu Chakrabarty. 2008. High-level synthesis of digital microfluidic biochips. *JETC* 3, 4 (2008). <https://doi.org/10.1145/1324177.1324178>
 - [50] Fei Su, William L. Hwang, and Krishnendu Chakrabarty. 2006. Droplet routing in the synthesis of digital microfluidic biochips, See [14], 323–328. <https://doi.org/10.1109/DATE.2006.244177>
 - [51] Ian I. Suni. 2008. Impedance methods for electrochemical sensors using nanomaterials. *TrAC Trends in Analytical Chemistry* 27, 7 (2008), 604 – 611. <https://doi.org/10.1016/j.trac.2008.03.012> Electroanalysis Based on Nanomaterials.
 - [52] Tao Xu and Krishnendu Chakrabarty. 2008. Integrated droplet routing and defect tolerance in the synthesis of digital microfluidic biochips. *JETC* 4, 3 (2008). <https://doi.org/10.1145/1389089.1389091>
 - [53] Tao Xu, Krishnendu Chakrabarty, and Fei Su. 2008. Defect-Aware High-Level Synthesis and Module Placement for Microfluidic Biochips. *IEEE Trans. Biomed. Circuits and Systems* 2, 1 (2008), 50–62. <https://doi.org/10.1109/TBCAS.2008.918283>

[54] Hailong Yao, Qin Wang, Yiren Shen, Tsung-Yi Ho, and Yici Cai. 2016. Integrated Functional and Washing Routing Optimization for Cross-Contamination Removal in Digital Microfluidic Biochips. *IEEE Trans. on CAD of Integrated Circuits and Systems* 35, 8 (2016), 1283–1296. <https://doi.org/10.1109/TCAD.2015.2504397>

[55] Ping-Hung Yuh, Chia-Lin Yang, and Yao-Wen Chang. 2007. Placement of defect-tolerant digital microfluidic biochips using the T-tree formulation. *JETC* 3, 3 (2007). <https://doi.org/10.1145/1295231.1295234>

[56] Ping-Hung Yuh, Chia-Lin Yang, and Yao-Wen Chang. 2008. BioRoute: A Network-Flow-Based Routing Algorithm for the Synthesis of Digital Microfluidic Biochips. *IEEE Trans. on CAD of Integrated Circuits and Systems* 27, 11 (2008), 1928–1941. <https://doi.org/10.1109/TCAD.2008.2006140>

[57] Yang Zhao and Krishnendu Chakrabarty. 2012. Cross-Contamination Avoidance for Droplet Routing in Digital Microfluidic Biochips. *IEEE Trans. on CAD of Integrated Circuits and Systems* 31, 6 (2012), 817–830. <https://doi.org/10.1109/TCAD.2012.2183369>